



# LMA ProSeal™ Proven versatility.




LMA ProSeal™ is a re-usable second generation, gastric access device which forms an effective First Seal™ with the oropharynx (oropharyngeal seal) and an innovative Second Seal™ with the upper oesophageal sphincter (oesophageal seal).

## The most versatile re-usable airway

- ▶ **Peace of mind** – Passive regurgitation can occur unexpectedly intraoperatively. LMA ProSeal™ enables the regurgitated fluid to pass up the drainage tube without leaking into the glottis<sup>1</sup>
- ▶ **Patient comfort** – LMA ProSeal™ reduces the likelihood of throat irritation and stimulation, and reduces post-operative nausea and vomiting by as much as 40% compared to an ETT<sup>2</sup>
- ▶ **Performance** – LMA ProSeal™ achieves a high seal pressure, with a median seal pressure of 32cm H<sub>2</sub>O<sup>3</sup>
- ▶ **Aspiration** – LMA ProSeal™ has a built-in drain tube that allows expelled gastric content to bypass the pharynx. This specific feature is designed to decrease the risk of aspiration<sup>4,5</sup>

Post anaesthesia care unit outcome data up to 24 hours<sup>2</sup>




★ Sore throat P<0.0001 vs. LMA ProSeal™ group.  
 ▲ Vomiting P<0.004 vs. LMA ProSeal™ group.  
 ♦ Nausea P<0.0001 vs. LMA ProSeal™ group.

Post-operative sore throat, nausea and vomiting were measured via patient interviews in a blind fashion.

“The LMA ProSeal™ currently has the broadest evidence to support its efficacy and safety profile.”

NAP4 report, 2011

The gastric drain tube enables confirmation of correct mask placement and sealing with the upper oesophageal sphincter. A patent drain tube reduces the risk of stomach insufflation.




**LMA ProSeal™: Product specification**

Mask size	Product code	Patient size	Maximum cuff volume (air)*	Largest size OG tube/salem pump
1	150010	Neonates/infants up to 5 kg	4 ml	2.7 mm / 8 Fr
1.5	150015	Infants 5-10 kg	7 ml	3.5 mm / 10 Fr
2	150020	Infants/children 10-20 kg	10 ml	3.5 mm / 10 Fr
2.5	150025	Children 20-30 kg	14 ml	4.9 mm / 14 Fr
3	150030	Children 30-50 kg	20 ml	5.5 mm / 16 Fr
4	150040	Adults 50-70 kg	30 ml	5.5 mm / 16 Fr
5	150050	Adults 70-100 kg	40 ml	6.0 mm / 18 Fr

\*These are maximum volumes that should never be exceeded. It is recommended that the cuff be inflated to a maximum of 60 cm H<sub>2</sub>O intracuff pressure. OG = orogastric

Distributed by:



PAB-2500-007 Rev B IW 201402

**TELEFLEX HEADQUARTERS INTERNATIONAL, IRELAND**  
 Teleflex Medical Europe Ltd., IDA Business and Technology Park,  
 Dublin Road, Athlone, Co Westmeath, Ireland  
 Phone +353 (0)9 06 46 08 00 Fax +353 (0)14 37 07 73  
 orders.intl@teleflex.com



Manufactured by:  
 The Laryngeal Mask Company Limited  
 Le Rocher, Victoria, Mahé, Seychelles



Consult IFU on this website:  
[www.LMACO.com](http://www.LMACO.com)

Teleflex [www.teleflex.com](http://www.teleflex.com) LMA [www.lmaco.com](http://www.lmaco.com)

Copyright © 2014 Teleflex Incorporated. All rights reserved. LMA, LMA ProSeal, First Seal, Second Seal, LMA Better by Design are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.

94 06 90 - AP 00 01 - REV B · LMA / PDF · 10 14 PDF

**References:**

1. Evans N.R. et al. CJA 2002; **49**: 413-416.
2. Hohlrieder M. et al. BJA 2007; **99**: 576-580.
3. Cook T.M., Gibbison B. BJA 2007; **99**: 436-439.
4. Mark D.A. CJA 2003; **50**: 78-80.
5. Brain A.I.J. et al. BJA 2000; **84**: 650-654.
6. 4th National Audit Project of the Royal College of Anaesthetists and the Difficult Airway Society: Major Complications of Airway Management in the United Kingdom. Report and findings: March 2011. Editors: Dr Tim Cook, Dr Nick Woodall and Dr Chris Frerk.



**Clinical evidence**



For the latest clinical evidence on LMA ProSeal™

[www.lmaco.com/evidence](http://www.lmaco.com/evidence)

**make-a-switch.com**



For more information on making the switch to second generation SADs

[www.make-a-switch.com](http://www.make-a-switch.com)



For the latest educational videos and clinician testimonials

[www.youtube.com/LaryngealMaskAirway](http://www.youtube.com/LaryngealMaskAirway)



For the latest news from LMA

[www.facebook.com/LMAInternational](http://www.facebook.com/LMAInternational)

**LMA™**  
 Better by Design



For product information and access to product instructions for use

[www.lmaco.com](http://www.lmaco.com)

**Teleflex**



For information on other products within the Teleflex product portfolio

[www.teleflex.com](http://www.teleflex.com)