


Pleur-evac™ Cactus™ and Atrium™ Oasis™ Series Dry Suction/Wet Seal Feature Comparison Overview

**Pleur-evac™ Cactus™
A-6000 Series**



**Atrium™ Oasis™
3600 Series**



SUCTION

Dry suction to -40 cmH₂O.

SUCTION SOURCE REQUIREMENTS

16 LPM

-80 mmHg

SUCTION CONTROL




Clicks in place at prescribed level.



Adjustable but does not click into place.

SUCTION INDICATOR



Distinct "On/Off" Orange Float only visible when suction has reached prescribed level (pre-set at -20 cmH₂O).

Variable bellows. Suction bellows will expand to the ▲ mark or beyond when suction is connected and operating at a regulator setting of -20 cmH₂O or higher. If the bellows is expanded but less than ▲ mark, increase the suction source vacuum to -80 mmHg or higher. For regulator settings less than -20 cmH₂O, and visible bellows expansion, the window will confirm suction operation.



SEAL

Water Seal must be established through the addition of 2 cm of water.
Patient is at risk of exposure to atmospheric air until water seal has been established.

PATIENT PROTECTION UPON TIP OVER


One-way valve protects the patient from atmospheric air in any orientation.
Suction ON: One-way valve protects the patient from atmospheric air in any orientation when suction is ON.
Suction OFF: Risk of atmospheric air entering the patient around angled fluid when suction is OFF.

Pleur-evac™ Cactus™ and Atrium™ Oasis™ Series

Dry Suction/Wet Seal Feature Comparison Overview



AIR LEAK METER



FLUID MIGRATION UPON TIP OVER

Fluid remains in seal when tipped. Drainage can migrate among sections of collection chamber.

COLLECTION CHAMBER

2500 ml

- 1 ml increments up to 100 ml, then 2 ml up to 200 ml
- Over 200 ml calibrations are in 5 ml increments

2100 ml

- 1 ml increments up to 100 ml, then 2 ml up to 200 ml
- Over 200 ml calibrations are in 10 ml increments

SPECIMEN SAMPLING

NEEDLELESS
Sample Port in patient tube.

NEEDLE
Not recommended

NEEDLE
20 gauge needle or smaller

INFANT UNIT

A-6020-08LF (1/2 ml increments up to 10 ml, then 1 ml to 150 ml)

3612 (1 ml up to 100 ml, then 2 ml up to 200 ml)

DUAL UNIT

A-6002-08LF (950 ml /1900 ml)

No – Discontinued

CONTINUOUS REINFUSION MODELS

A-9250LF (2100 ml; first 1000 ml recoverable for reinfusion)

No – Discontinued

Atrium and Oasis are trademarks of Atrium Medical Corporation. Federal Law (USA) restricts these devices to sale by or on the order of a physician. Not all products may be available in all countries. For product information please contact la.cs@teleflex.com or your local representative. Refer to the applicable Instructions for Use for the indications approved in your geography. Please check your local regulatory approval status. Teleflex, the Teleflex logo, Cactus, and Pleur-evac are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Other names may be trademarks of their respective owners. Revised: 09/2023. © 2023 Teleflex Incorporated. All rights reserved. MC-009246 LA EN