


3-WAY SILICONE FOLEY CATHETER FOR POSTOPERATIVE URINARY DRAINAGE

Silicone Foley catheter for optimum patient care

Our well-known Rüsçh 3-way silicone Foley catheters are first and foremost distinguished by their smooth surface, both on the outer and inner lumen, as well as by their balloon which fits closely to the shaft. Urethral irritation and encrustation on the catheter inside are thus largely avoided.



Because of the material's good tolerability, these silicone Foley catheters are especially suited for long-term urinary drainage, and can also be used with extremely sensitive patients.

Rüsçh 3-way silicone Foley catheters offer sufficient stability for simple and safe application, and the spacious inner lumen allows high drainage and irrigation capacity. Yet, they remain flexible and soft to the touch for maximum patient comfort.



AT A GLANCE





- transparent 3-way silicone Foley catheters for postoperative drainage
- closely fitting balloon with large volume
- colour-coded valve for Luer and Luer-lock syringe tips
- white radiopaque tip and white contrast stripe
- catheter length approx. 42 cm
- sterile, single use
- latex-free

3-WAY SILICONE FOLEY CATHETER FOR POSTOPERATIVE URINARY DRAINAGE





SILICONE FOLEY CATHETER, TRANSPARENT

RUSCH

SIZE	BALLOON	QTY	REF.	CATHETER FOR POSTOPERATIVE DRAINAGE 3-WAY
 Cylindrical, 2 staggered eyes				<ul style="list-style-type: none"> • valve for Luer and Luer-lock syringe tips • white radiopaque tip and white contrast stripe • length approx. 42 cm • sterile, single use • latex-free
Ch. 16	40 ml	5	173800	
Ch. 18	50 ml	5		
Ch. 20	60 ml	5		
Ch. 22	70 ml	5		
Ch. 24	80 ml	5		
 Couvelaire, 2 staggered eyes				
Ch. 16	40 ml	5	173900	
Ch. 18	50 ml	5		
Ch. 20	60 ml	5		
Ch. 22	70 ml	5		
Ch. 24	80 ml	5		
 Dufour, 2 staggered eyes				
Ch. 16	40 ml	5	174000	
Ch. 18	50 ml	5		
Ch. 20	60 ml	5		
Ch. 22	70 ml	5		
Ch. 24	80 ml	5		
 Mercier, 2 staggered eyes				
Ch. 16	40 ml	5	174100	
Ch. 18	50 ml	5		
Ch. 20	60 ml	5		
Ch. 22	70 ml	5		
Ch. 24	80 ml	5		

SIZES Ch. 16 Ch. 18 Ch. 20 Ch. 22 Ch. 24



Distributed by:
TELEFLEX MEDICAL Europe • IDA Business Park • Athlone • Co Westmeath • Ireland
 Tel. +353 (0)9 06 46 08 00 • Fax +353 (0)14 37 07 73 • orders.intl@teleflexmedical.com
Teleflex Medical UK • High Wycombe • Bucks. HP 12 3 ST • United Kingdom
 Tel. +44 14 94 53 27 61 • Fax +44 14 94 52 46 50 • info.uk@teleflexmedical.com
 www.ruesch.de • www.teleflexmedical.com