


**The Suture of Choice—
Deklene® II**

- Minimal memory provides superior handling, accurate knot placement, and knot security
- Uniform diameter with high tensile strength resists breakage
- Soft, supple material reduces memory and improves handling
- Available in a variety of packaging configurations, including straight packs, packs of one, two and four sutures
- All packaging is conveniently designed for direct-arm access



Convenient direct-arm access. Large selection of needle/suture combinations is available.



the suture of choice—
Deklene® II

Ordering Information

A wide assortment of Deklene® II suture/needle combinations is available for cardiovascular and other specialty surgeries. For more information, contact your Teleflex Medical Surgical Specialist or call 800-367-7874.




Teleflex Medical
P.O. Box 12600
2917 Weck Drive
Research Triangle Park, NC 27709

USA Phone: (800) 367-7874
USA Fax: (800) 677-6688

International Phone: (919) 361-3961
International Fax: (919) 361-4111
E-Mail: international@teleflexmedical.com

www.teleflex.com

Deklene® II *Polypropylene Suture* XT Series *Premium Needles*




The suture of choice combined with the needle of choice for coronary artery bypass grafting and other specialty surgical procedures.

the needle of choice— *Premium XT Series*

At the moment the bypass graft is sutured in place, exceptional needle performance is critical. Suture handling and security become essential. Teleflex Medical understands that the perfect CABG needle and suture combination is a natural extension of skilled surgical hands.




The Needle of Choice— Premium XT Series




- 300 series steel blends superior strength and ductility
- Extraordinary point geometries enable easy penetration of calcified vessels
- Square-body design enhances strength and offers a closer needle/suture ratio
- Strong laser-drilled attachment increases needle area available to grasp with needle holders without the risk of bending
- Available on Deklene® II sizes 5/0 through 8/0






Extraordinary point geometries and
custom 300 series steel.




NEEDLE/SUTURE COMBINATIONS


Needle	No. of Needles	Suture	Length x No. of Sutures	Code No.
XT-7 	2	8/0	18"	D-7501A
	2	8/0	18"x1	D-7501M1A
	2	8/0	18"x2	D-7501M2
	2	8/0	18"x4	D-7501M4
 <i>Taper Point</i>	2	8/0	24"	D-7504
	2	8/0	24"x1	D-7504M1
	2	8/0	24"x2	D-7504M2
	2	8/0	24"x4	D-7504M4
	2	8/0	30"x1	D-6551M1



XT-6.5 	2	8/0	24"	D-7505A
	2	8/0	24"x1	D-7505M1A
	2	8/0	24"x2	D-7505M2
	2	8/0	24"x4	D-7505M4
 <i>Taper Point</i>	2	7/0	24"	D-7506A
	2	7/0	24"x1	D-7506M1
	2	7/0	24"x2	D-7506M2
	2	7/0	24"x4	D-7506M4
	2	7/0	30"x2	D-7507M2
	2	7/0	30"x4	D-7507M4


XT-6 	2	7/0	18"	D-7062A
	2	7/0	18"x4	D-7062M4A
	2	7/0	24"	D-7004A
	2	7/0	24"x1	D-7004M1A
 <i>Taper Point</i>	2	7/0	24"x4	D-7004M4A
	2	7/0	30"x2	D-7063M2A
	2	7/0	30"x4	D-7063M4A
	2	6/0	18"	D-793A
	2	6/0	18"x4	D-793M4A
	2	6/0	24"	D-5007A
	2	6/0	24"x4	D-5007M4A
	2	6/0	30"	D-7070A
2	6/0	30"x4	D-7070M4A	

XT-6 	1	7/0	18"	D-790A
 <i>Taper Point</i>				

XT-1 	2	6/0	18"	D-655A
	2	6/0	18"x4	D-655M4A
	2	6/0	24"	D-7016A
	2	6/0	24"x2	D-7016M2A
 <i>Taper Point</i>	2	6/0	24"x4	D-7016M4A
	2	6/0	30"	D-7076A
	2	6/0	30"x1	D-7076M1A
	2	6/0	30"x2	D-7076M2A
	2	6/0	30"x4	D-7076M4A
	2	6/0	36"	D-7017A
	2	6/0	36"x4	D-7017M4A
	2	5/0	24"	D-7060A



Needle	No. of Needles	Suture	Length x No. of Sutures	Code No.
XT-1 <i>continued</i> 	2	5/0	24"x4	D-7060M4A
	2	5/0	30"	D-7018A
	2	5/0	36"	D-7061A
	2	5/0	36"x2	D-7061M2A
	2	5/0	36"x4	D-7061M4A


XT-1 	1	6/0	30"	D-730A
 <i>Taper Point</i>				




EZ-1 	2	6/0	30"	D-7414
	2	6/0	30"x2	D-7414M2
	2	6/0	30"x4	D-7414M4

 <i>Taper Cutting</i>				
--	--	--	--	--

EZ-6 	2	7/0	24"x2	D-7405M2
				
 <i>Taper Cutting</i>				

T-1 	2	7/0	24"x4	D-7015M4
	2	5/0	18"	D-660
				

 <i>Taper Point</i>				
--	--	--	--	--

ST-6 	2	8/0	24"x2	D-7308M2
				
 <i>Slim Taper Point</i>				