


DEKNATEL | WECK

Dekna-lok is designed to facilitate soft tissue approximation during open and laparoscopic surgery.

Dekna-lok is comprised of two non-absorbable polymer components affixed to a Size 1 absorbable polyglycolic acid suture with a 39mm taper needle. The polymer components are designed to spread suture knot pressure across tissue.



ORDERING INFORMATION

DESCRIPTION	ITEM NO.
Dekna-lok	BP1000V2L

CONTACT INFORMATION

United States
 Phone: (866) 246-6990 Fax: (800) 932-5329

Europe
 Phone: 33 562 18 79 10 Fax: 33 561 83 35 84

international@teleflexmedical.com

Some products may not be available in all markets. Please contact your Teleflex Medical distributor for further details.

For mailing addresses, please visit our web page at: www.teleflexmedical.com/contact.html

Weck® is a registered trademark of Teleflex Inc. Dekna-lok™ is a trademark and Deknatel® is a registered trademark of Teleflex-CT Devices Incorporated.

© 2006 Teleflex Medical All rights reserved. Lit. No. 25-4002-EN 2006-06



Teleflex Medical 2917 Weck Drive
 Research Triangle Park, NC 27709
 USA Tel: 866-246-6990
 Int'l: 919-361-3961
www.teleflexmedical.com

