

NEW

ARROW™ EZ-IO™
Intraosseous Access Procedure Tray



Arrow™ EZ-IO™ Intraosseous Access Procedure Tray

In urgent and emergent situations, every moment matters. The Arrow™ EZ-IO™ Intraosseous Access Procedure Tray has everything healthcare professionals need to gain intraosseous access at the speed of life.

Complete, All-In-One Tray

- Helps streamline workflow and increase efficiency
- Single-use and sterile, appropriate for multiple care settings
- Supports standardization across the hospital

Clinical Convenience

- Compact design for storage on hospital shelving
- Small enough to fit in any code cart for increased accessibility
- Includes a step-by-step Quick Reference Guide

Industry Tested & Trusted

- 3 seconds for fluid and medication delivery to the heart through the proximal humerus^{1*}
- 97% first attempt success rate²
- Less than 1% serious complication rate³

Benefits



The Clinician

Convenient tray design provides IO procedure components in one place



The Institution

Sterile, pre-packaged tray supports standardization and workflow efficiency



The Patient

IO access availability for all difficult vascular access patients across multiple care settings

Teleflex™

Get to know the EZ-IO™ Procedure Tray:




EZ-IO™ Single-Use Power Driver

Sterile, battery-powered driver for convenient IO insertion



EZ-IO™ Needle Set

Includes one 15, 25, or 45 mm needle, proven to deliver fast,⁴ precise, and steady insertion



... ChloroPrep™ One-Step Frepp™ Solution

Broad spectrum antiseptic for prepping the skin



... Arrow SharpsAway™ II Locking Disposal Cup

Promotes sharps safety




... EZ-Connect™ Extension Set

Reinforced 90° angle to help prevent your line from kinking



... Pre-filled Saline Syringe

To perform a rapid flush to displace marrow and facilitate flow



... EZ-Stabilizer™ Dressing

Secures placement and protects the insertion site



... EZ-IO™ Wristband

To document insertion time

EZ-IO™ Procedure Tray Ordering Information

DESCRIPTION	ITEM NUMBER	PATIENT WEIGHT	QTY/CASE
15 mm Arrow™ EZ-IO™ Intraosseous Access Tray + Single-Use Power Driver	9018T-VC-005	3-39 kg	5
25 mm Arrow™ EZ-IO™ Intraosseous Access Tray + Single-Use Power Driver	9001T-VC-005	≥3 kg	5
45 mm Arrow™ EZ-IO™ Intraosseous Access Tray + Single-Use Power Driver	9079T-VC-005	≥40 kg	5

*Based on adult proximal humerus study conducted in a healthy individual.

References:

1. Montez DF, Puga T, Miller L, et al. *Ann Emerg Med.* 2015;66(4S):S47. Research sponsored by Teleflex Incorporated.
2. Brenner T, Bernhard M, Helm M, et al. *Resuscitation.* 2008;78(3):314-319.
3. Teleflex Internal Data on File 2018.
4. Davidoff J, Fowler R, Gordon D, et al. *JEMS.* 2005;30(10):s20-s23. Research sponsored by Teleflex Incorporated.

The Arrow™ EZ-IO™ Intraosseous Access Procedure Tray is indicated for intraosseous access anytime in which vascular access is difficult to obtain in emergent, urgent or medically necessary cases for up to 24 hours. For patients ≥ 12 years old, the device may be extended for up to 48 hours when alternate intravenous access is not available or reliably established.

Rx Only.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

The Arrow™ EZ-IO™ Needle Set is Sterile, Single Use: Do not reuse, reprocess or re-sterilize. Reuse of device creates a potential risk of serious injury and/or infection, which may lead to death. Refer to Instructions for Use for complete warnings, indications, contraindications, precautions, and potential complications.

ChloroPrep and Frepp are trademarks or registered trademarks of Carefusion or one of its subsidiaries.

Teleflex, the Teleflex logo, Arrow, EZ-IO, EZ-Connect, EZ-Stabilizer, and SharpsAway are trademarks or registered trademarks of Teleflex Incorporated and its affiliates, in the U.S. and/or other countries. All other trademarks or registered trademarks are the property of their respective owners. All rights reserved. ©2024 Teleflex Incorporated. MC-009484 (01/2024)