

RÜSCH

HANS RUDOLPH V2
REUSABLE
VENTILATION MASKS



Teleflex
MEDICAL

REUSABLE VENTILATION MASKS HANS RUDOLPH V2

Hans Rudolph Version 2 Oronasal Mask Non-invasive ventilation (NIV)

- Transparent mask made of 100% silicone
- Anatomically shaped silicone discs facilitate flexible positioning
- Thin silicone membrane in the nasal area minimizes leaks
- Shaped chin section for secure fit
- Smallest possible dead space
- Latex-free
- Manufactured from a single piece
- Connector permits 360° rotation
- Reusable




When ordering,
please always
indicate
Ref. and size.

REF.	DESIGNATION	DESCRIPTION	SIZES	QTY
15 43 17	V2 mask (NIV)	Oronasal mask with green elbow connector 22 mm I.D. and headgear	P/ES/S/M/L	1
15 43 68	V2 mask	Mask only, without connector	P/ES/S/M/L	1
15 43 25	Headgear, large (L)	suitable for mask size L		1
15 43 25	Headgear, medium (M)	suitable for mask sizes M and S		1
15 43 25	Headgear, small (S)	suitable for mask sizes P and ES, alternatively also for larger masks		1
15 43 73	Green elbow connector	22 mm I.D. (NIV)		1
15 43 74	Green connector	straight, 22 mm I.D. (NIV)		1
15 43 69	Replacement clips	for headgear		4
15 43 38	Silicone sealing caps	for connector		2
15 43 78	Mask sizing gauge	to determine correct mask size		1

Metro Seal

REF.	DESIGNATION	QTY
15 43 70	Metro Seal	4

- To lead through probes and cables (Ch. 4 - 18)
- Additional sealing and adaptation for difficult anatomical conditions
- Made from biotolerant transparent hydrogel



Secure hold and high patient comfort

Hans Rudolph Version 2 Oronasal Mask CPAP / BiLevel

- Blue elbow connector with anti-asphyxia valve (AAV); the safety valve opens at a pressure below 3 cm H₂O
- With exhalation opening in the connector



REF.	DESIGNATION	DESCRIPTION	SIZES	QTY
15 43 22	V2 mask (CPAP)	Oronasal mask with blue elbow connector 22 mm O.D., with AAV and exhalation opening, headgear	P/ES/S/M/L	1
15 43 68	V2 mask	Mask only, without connector	P/ES/S/M/L	1
15 43 25	Headgear, large (L)	Suitable for mask size L		1
15 43 25	Headgear, medium (M)	Suitable for mask sizes M and S		1
15 43 25	Headgear, small (S)	Suitable for mask sizes P and ES, alternatively also for larger masks		1
15 43 77	Blue elbow connector	22 mm O.D. (CPAP) with AAV and exhalation opening		1
15 43 69	Replacement clips	For headgear		4
15 43 38	Silicone sealing caps	For connector		2
15 43 78	Mask sizing gauge	To determine correct mask size		1

*When ordering,
please always
indicate
Ref. and size.*

Hans Rudolph Version 2 Oronasal Mask BiLevel

- Yellow elbow connector with anti-asphyxia valve (AAV); the safety valve opens at a pressure below 3 cm H₂O
- Without exhalation opening in the connector



REF.	DESIGNATION	DESCRIPTION	SIZES	QTY
15 43 20	V2 mask (BiLevel)	Oronasal mask with yellow elbow connector 22 mm I.D., with AAV and headgear	P/ES/S/M/L	1
15 43 68	V2 mask	Mask only, without connector	P/ES/S/M/L	1
15 43 25	Headgear, large (L)	Suitable for mask size L		1
15 43 25	Headgear, medium (M)	Suitable for mask sizes M and S		1
15 43 25	Headgear, small (S)	Suitable for mask sizes P and ES, alternatively also for larger masks		1
15 43 79	Yellow elbow connector	22 mm I.D. (BiLevel) mit AAV, ohne Ausatemöffnung		1
15 43 69	Replacement clips	For headband		4
15 43 38	Silicone sealing caps	For connector		2
15 43 78	Mask sizing gauge	To determine correct mask size		1

*When ordering,
please always
indicate
Ref. and size.*

ARROW

RÜSCH

 HUDSON RCI

 Gibeck

RÜSCH CARE

MEDICAL  SERVICE

WECK

Pilling

PLEUR-EVAC



TAUT

Teleflex serves a wide variety of customers around the world by designing, manufacturing, and distributing specialty-engineered products. Distinguished by a significant presence in healthcare, Teleflex also serves the aerospace and industrial markets with businesses that have significant market position and brand recognition.

Teleflex Medical is a global leader in disposable medical products for critical care and surgical applications and a business segment of Teleflex Incorporated (NYSE: TFX). The company also produces surgical instruments and devices, cardiac devices and other specialty products for device manufacturers. The Teleflex Medical family of brands includes the trusted names of Arrow®, Beere®, Deknatel®, Gibeck®, HudsonRCI®, KMedic®, Pilling®, Pleur-Evac®, Rüschi®, Sheridan®, SMD, SSI™, Taut™, TFX Medical OEM and Weck®.

www.ruesch.de
www.teleflex.com

Manufacturer:
Hans Rudolph, inc. 8325, Cole Parkway,
Shawnee, Kansas, 66227, USA
ISO 9001 | ISO 16485 | EN 46001

Distributed by:
NETHERLANDS TELEFLEX MEDICAL B.V., P.O. Box 40, 3640 AA Mijdrecht,
Phone +31 (0)297 540 703, Fax +31 (0)297 533 726, info.nl@teleflexmedical.com
UNITED KINGDOM TELEFLEX MEDICAL, High Wycombe, Bucks. HP12 3ST
Phone +44 (0)1494 53 27 61, Fax +44 (0)1494 52 46 50, info.uk@teleflexmedical.com