

## HUDSON RCI SOFTECH® PLUS NASAL CANNULA

Soft. Safe. Superior.



# HUDSON RCI SOFTECH® PLUS NASAL CANNULA

## BREATHING, MADE EASIER®

Trusted brands make Teleflex a reliable and strong partner. Built on a solid tradition of innovation, Teleflex is a global leader in superior medical supplies designed to help providers minimise risk and maximise outcomes for their patients.

When it comes to respiratory care, patients and caregivers around the world look to HUDSON RCI for innovative products that provide added safety and comfort. Guided by the voice of the respiratory clinician, we strive not only to advance patient care, but also to raise caregiver efficiency.

As a result of pressure, the use of cannulas can cause discomfort among patients. To remove some of the pressure, we are pleased to introduce the DEHP-free Softech Plus nasal cannula from HUDSON RCI.

## TELEFLEX – HIGH QUALITY MEDICAL SUPPLIES FROM A SINGLE SOURCE

## OUR SOFTEST NASAL CANNULA EVER

Softech Plus is the ideal choice for comfortable oxygen delivery because of the combination of a softer DEHP-free material and improved product design. The use of Softech Plus will help improve patient comfort and enhance the quality of care.

### SOFT

Thanks to a special new blend of materials, Softech Plus is our softest nasal cannula ever, and ideally fits most patients, from adults to neonates. Extending from the nasal prongs to the entire lariat tubing, the super soft material can help reduce surface friction by 30% compared to other cannulas\*. Featuring an over-the-ear-design, this cannula ensures optimal patient comfort for the nose, the face and the ears.


### SAFE

In fact, the Softech Plus nasal cannula offers a new level of patient safety. The use of **DEHP-free materials** helps clinicians reduce their patients' exposure to DEHP, while the **kink-resistant star-lumen lariat tubing** (adult and pediatric sizes only) ensures the uninterrupted flow of oxygen to the patient. In addition, the highly visible, green-coloured oxygen tubing may prevent connection errors and accidental tripping.

### SUPERIOR

Last but not least, we have integrated a special notch between the nasal prongs (except item code 1873) to achieve greater flexibility in terms of a more individual fit. Moreover, Teleflex provides a comprehensive range of sizes, including adult, paediatric, infant and neonatal sizes.





*\*Friction can contribute to pressure ulcer formation.<sup>1</sup> In laboratory tests, the HUDSON RCI Softech Plus cannula demonstrated 30% less surface friction compared to alternative cannulas in its category.<sup>2</sup>*

**References**

1 International review. Pressure ulcer prevention: pressure, shear, friction and microclimate in context. A consensus document.

London: Wounds International, 2010.

2 Data on file with Teleflex Incorporated.

**KEY BENEFITS**

Optimal patient comfort:

- reduced surface friction
- flexible fit
- comprehensive range includes cannulas for neonates

Increased safety:

- uninterrupted oxygen flow (kink-resistant star lumen lariat in adult and paediatric sizes)
- reduced exposure to DEHP
- may reduce connection errors and accidental tripping

**SOFTECH PLUS NASAL CANNULA**

**HUDSON RCI**

REF.	SIZE	TUBING LENGTH (CM)	STAR LUMEN® TUBING	STANDARD CONNECTOR	UNIVERSAL CONNECTOR	QTY
1852	adult	420	•		•	50
1853	paediatric	420	•		•	50
1870	adult	210	•	•		50
1871	paediatric	210	•	•		50
1872	infant	210	•	•		50
1873	neonatal	210	•	•		50
1874	adult	210	•		•	50
1876	paediatric	210	•		•	50
1877	adult	420	•	•		50
1878	paediatric	420	•	•		50

Teleflex is a leading global provider of specialty medical devices used for diagnostic and therapeutic procedures in critical care, urology and surgery. Our mission is to provide solutions that enable healthcare providers to improve outcomes and enhance patient and provider safety. We specialise in devices for general and regional anaesthesia, cardiac care, respiratory care, urology, vascular access and surgery and we serve healthcare providers in more than 150 countries. Teleflex also provides specialty products for medical device manufacturers.

Our well known brands include ARROW®, DEKNATEL®, GIBECK®, HUDSON RCI®, KMEDIC®, LMA™, PILLING®, PLEUR-EVAC®, RÜSCH®, SHERIDAN®, TAUT®, TFX OEM®, VASONOVA™, VIDACARE® and WECK®, all of which are trademarks or registered trademarks of Teleflex Incorporated.

Teleflex global operations: Australia, Austria, Belgium, Canada, China, Czech Republic, France, Germany, Greece, India, Ireland, Italy, Japan, Malaysia, Mexico, Netherlands, New Zealand, Portugal, Singapore, Slovak Republic, South Africa, Spain, Switzerland, United Kingdom, Uruguay and USA.

#### **YOUR INTERNATIONAL CONTACTS:**

##### **TELEFLEX HEADQUARTERS INTERNATIONAL, IRELAND**

Teleflex Medical Europe Ltd., IDA Business and Technology Park,  
Dublin Road, Athlone, Co Westmeath  
Phone +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73  
orders.intl@teleflex.com

**AUSTRALIA/NEW ZEALAND** +61 (0)3 9081 0600

**AUSTRIA** +43 (0)1 402 47 72

**BELGIUM** +32 (0)2 333 24 60

**CHINA (SHANGHAI)** +86 (0)21 6163 0965

**CHINA (BEIJING)** +86 (0)10 6418 5699

**CZECH REPUBLIC** +420 (0)495 759 111

**FRANCE** +33 (0)5 62 18 79 40

**GERMANY** +49 (0)7151 406 0

**GREECE** +30 210 67 77 717

**INDIA** +91 (0)44-2836 5040

**ITALY** +39 0362 58 911

**JAPAN** +81 (0)3 3379 1511

**NETHERLANDS** +31 (0)88 00 215 00

**PORTUGAL** +351 22 541 90 85

**SINGAPORE** +65 6439 3000

**SLOVAK REPUBLIC** +421 (0)3377 254 28

**SOUTH AFRICA** +27 (0)11 807 4887

**SPAIN** +34 918 300 451

**SWITZERLAND** +41 (0)31 818 40 90

**UNITED KINGDOM** +44 (0)1494 53 27 61

For detailed information see [www.teleflex.com](http://www.teleflex.com)

The products in this catalogue may not be available in all countries. Please contact your local representative. All data current at time of printing (09/2014). Subject to technical changes without further notice.

94 06 73 - 00 00 01 · REV C · MC / WM · 09 14 01