



RÜSCH CRYSTALCLEAR PLUS

Better comfort and more safety for user and patient



AT A GLANCE

- single-use system for infection prevention
- flexible inner cannula system for ease of handling
- colour-coded tube neckplate and inner cannula stopper for quick and safe size identification
- pilot balloon valve sealing cap
- thermosensitive tube material for greater flexibility
- soft, open, transparent fixation flange for comfortable fit
- continuous X-ray marker along the tube shaft
- integrated rotating 15 mm connector
- sterile, latex-free, for single use

COMFORT & EFFICIENCY

Our CrystalClear Plus tracheostomy tube set is an advanced tube system with a flexible inner cannula system. It is distinguished by top quality, cost effectiveness and greatest wearing comfort. Simple and safe to use it renders the care of tracheostomy patients comfortable and efficient.

SAFETY & EASE OF HANDLING

CrystalClear Plus is a single-use system which helps to reduce the risk of infection to a minimum, as the specially designed inner cannula system protects the inner tube lumen from secretion. Removal of secretions is simple and efficient at an early stage of accumulation through frequent changes of the inner cannula,

and well-aligned accessory parts assist the patient in clearing mucus. With cuffed tubes, a pilot balloon valve cap is premounted to protect the valve and the inflated cuff.

BROAD SIZE RANGE

Our CrystalClear Plus is available from sizes I.D. 3.5 to 10.0 mm, in 0.5 mm steps. The optimum I.D. and O.D. ratio between tube and inner cannula ensures best possible results, both in paediatric as well as adult applications.



CRYSTALCLEAR PLUS SET (with cuff, Ref. 121303, without cuff, Ref. 121300)
RUSCH

	ORDER SIZE / I.D. IN MM	O.D. IN IN MM	CLL* IN MM	DIM. A IN MM	DIM. B IN MM	DIM. C IN MM	ANGLE θ APPROX.	CUFF Ø IN MM **	INNER CANNULA	QTY
Paediatrics	● 3.5	5.3	47	17	0	30	95°	11	–	5
	○ 4.0	5.9	54	18	0	36	95°	12	–	5
	● 4.5	6.5	54	19	0	35	95°	14	–	5
	● 5.0	7.0	57	20	0	37	95°	14	–	5
	● 5.5	7.5	58	20	0	38	95°	15	–	5
	● 6.0	8.7	64	25.5	0	38.5	95°	19	6.0	5
Adult	● 6.5	8.7	64	25.5	0	38.5	95°	19	6.5	5
	● 7.0	10.0	70	29.2	0	40.8	95°	21	7.0	5
	● 7.5	10.7	75.5	30	0	45.5	95°	25	7.5	5
	● 8.0	11.3	81	35.5	0	45.5	95°	28	8.0	5
	● 8.5	11.3	81	35.5	0	45.5	95°	28	8.5	5
	● 9.0	12.7	87	39.9	0	47.1	95°	28	9.0	5
	● 9.5	12.7	87	39.9	0	47.1	95°	28	9.5	5
	● 10.0	14.0	95.5	48.4	0	47.1	95°	30	10.0	5

* Total length (CLL=A+B+C)

** only 121303 CrystalClear Plus Set

CRYSTALCLEAR PLUS PDT SET (Ref. 121304)
RUSCH




ORDER SIZE / I.D. IN MM	O.D. IN IN MM	CLL* IN MM	DIM. A IN MM	DIM. B IN MM	DIM. C IN MM	ANGLE θ (APPROX.)	CUFF Ø IN MM	QTY
● 7.0	10.0	70.0	29.2	–	40.8	95°	21	1
● 8.0	11.3	81.0	35.5	–	45.5	95°	28	1
● 9.0	12.7	87.0	39.9	–	47.1	95°	28	1

* Total length (CLL=A+B+C)

**CRYSTALCLEAR PLUS
INNER CANNULA (Ref. 121310)**
RUSCH

ORDERSIZE	ID IN MM	OD IN MM	QTY
● 6.0	5.0	5.7	10
● 6.5	5.0	5.7	10
● 7.0	6.0	6.7	10
● 7.5	6.5	7.3	10
● 8.0	7.0	7.7	10
● 8.5	7.0	7.7	10
● 9.0	8.0	8.7	10
● 9.5	8.0	8.7	10
● 10.0	8.8	9.7	10

ACCESSORIES
RUSCH

	REF.	DESCRIPTION	QTY
	121311	shower cap for CrystalClear Plus	4
	121312	cough cap for CrystalClear Plus	4
	121313	sealing cap for CrystalClear Plus	4

**CRYSTALCLEAR PLUS
SET COMPONENTS**

- CrystalClear Plus tracheostomy tube with or without low-pressure cuff and premounted valve sealing cap
- 2 flexible inner cannulas with colour-coded stopper, starting from size 6
- insertion aid (obturator), starting from size 6
- soft tracheostomy tube neckband
- water-soluble lubricating jelly, 3g sachet

**CRYSTALCLEAR PLUS PDT
SET COMPONENTS**

- tracheostomy tube with low-pressure cuff
- flexible fixation flange
- comfort neck band
- introduction dilator
- sterile
- single use

DISTRIBUTED BY:

TELEFLEX HEADQUARTER INTERNATIONAL, IRELAND · Teleflex Medical Europe Ltd. · IDA Business & Technology Park
Dublin Road · Athlone · Co Westmeath · Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com

UNITED KINGDOM Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com

WWW.TELEFLEX.COM

94 05 59 - 00 00 01 · REV B · 09 13 01 · All data current at time of printing. Subject to technical changes without further notice.

teleflex