



## RÜSCH SAFETY SILK SERIES

A touch of silk



**Teleflex**


# A TOUCH OF SILK


For your most demanding patients, as well as for everyday use, Teleflex has designed a new range of endotracheal tubes. Extremely soft and flexible yet kink-resistant, with an optimal ID/OD ratio, the Safety Silk Series' material adapts to every patient's anatomy. It is particularly well suited for nasal intubation where softness counts and for more delicate paediatric patients. When your patient needs extra care, our Safety Silk Series provides for a smooth, atraumatic and easy intubation.

Your concerns are ours. Most of our Safety Silk products have been re-engineered to bring you the first ever available DEHP-free endotracheal tubes on the market. By increasing overall safety, we trust this new development provides an answer to some of your concerns for your patients' well-being.

But not only with regard to DEHP-free endotracheal tubes have we listened to you. We have also extended the proximal length of our Nasal Safety Silk products to avoid trauma to patients' foreheads while fitting. This new length together with all the well-known advantages of our Nasal Safety Silk range (softness, comfort, optimal curvature) are additional reasons why you should choose the Teleflex brand.

**SAFETY SILK SERIES BY TELEFLEX:  
A TOUCH OF SILK FOR ALL YOUR PATIENTS.**

 ALL TUBES IN THE SAFETY SILK RANGE ARE AVAILABLE IN PAEDIATRIC SIZES



## SUPER SAFETY SILK

CLEAR ET TUBE MADE OF SOFT DEHP-FREE PVC WITH CUFF

- nasal/oral
- 1 eye (Murphy)
- tissue-friendly thermosensitive material, especially advantageous for long-term intubation
- Safety Silk surface
- kink-resistant
- soft atraumatic tip
- low pressure cuff
- bonded connector
- continuous blue X-ray marker
- graduated
- blue pilot balloon
- inflation channel exit close to connector
- black ring for intubation depth control
- valve for Luer and Luer-lock syringes
- latex-free
- sterile
- single use




Super Safety Silk ensures easy and safe insertion (adult nasal & paediatric) and is perfectly suited for long-term intubation.

## NASAL SAFETY SILK

CLEAR ANATOMICALLY SHAPED ET TUBE MADE OF SOFT DEHP-FREE PVC WITH CUFF

- nasal
- 1 eye (Murphy)
- tissue-friendly thermosensitive material, especially advantageous for long-term intubation
- Safety Silk surface
- kink-resistant
- soft atraumatic tip
- low-pressure cuff
- bonded connector
- continuous blue X-ray marker
- graduated
- blue pilot balloon
- inflation channel exit close to curve to avoid patient trauma when tube trimming is required
- black ring for intubation depth control
- valve for Luer and Luer-lock syringes
- latex-free
- sterile
- single use




Nasal Safety Silk is specifically designed for atraumatic insertion and smoothly adapts to suit the nasal anatomy.

## PEDI SAFETY SILK

CLEAR ET TUBE MADE OF SOFT PVC  
WITHOUT CUFF

- nasal/oral (Magill)
- extremely soft body & tip material for atraumatic insertion in delicate patients
- tissue-friendly thermosensitive material, especially advantageous for long-term intubation
- Safety Silk surface
- kink-resistant
- black tip for easy positioning
- semi-seated connector
- continuous blue X-ray marker
- graduated
- latex-free
- sterile
- single use



Pedi Safety Silk is the tube of choice for uncuffed intubation in children & neonates.

## NASOPHARYNGEAL SAFETY SILK

CLEAR NASOPHARYNGEAL AIRWAY MADE OF SOFT DEHP-FREE PVC

- tissue-friendly, thermosensitive material
- Safety Silk surface
- kink-resistant
- soft atraumatic tip
- latex-free
- sterile
- single use



The nasopharyngeal Safety Silk is the NPA of choice for an easy & atraumatic insertion, providing airway patency in all situations.

ALL TUBES IN THE SAFETY SILK SERIES ARE AVAILABLE IN PAEDIATRIC SIZES.

**SUPER SAFETY SILK**

**RÜSCH**

REF.		ORDER SIZE/I.D.	O.D.	CUFF Ø	LENGTH	QTY
112682 - Murphy eye	Paediatric	3.0 mm	5.0 mm	8.0 mm	170 mm	10
		3.5 mm	5.3 mm	10.5 mm	190 mm	10
		4.0 mm	6.0 mm	10.5 mm	220 mm	10
		4.5 mm	6.3 mm	13.0 mm	230 mm	10
112680 - Magill 112682 - Murphy eye	Adult	5.0 mm	6.7 mm	16.5 mm	250 mm	10
		5.5 mm	7.3 mm	18.5 mm	280 mm	10
		6.0 mm	8.0 mm	20.5 mm	290 mm	10
		6.5 mm	8.7 mm	24.0 mm	300 mm	10
		7.0 mm	9.3 mm	26.0 mm	320 mm	10
		7.5 mm	10.0 mm	26.0 mm	330 mm	10
		8.0 mm	10.7 mm	28.0 mm	340 mm	10
		8.5 mm	11.3 mm	28.0 mm	345 mm	10
		9.0 mm	12.0 mm	28.0 mm	350 mm	10
		9.5 mm	12.7 mm	29.0 mm	350 mm	10
		10.0 mm	13.3 mm	29.0 mm	350 mm	10



**NASAL SAFETY SILK**

**RÜSCH**

REF.		ORDER SIZE/I.D.	O.D.	CUFF Ø	TOTAL LENGTH	QTY
111782	Paediatric	3.0 mm	5.3 mm	11.0 mm	245 mm	10
		3.5 mm	6.0 mm	11.0 mm	271 mm	10
		4.0 mm	6.3 mm	13.0 mm	303 mm	10
		4.5 mm	6.7 mm	13.0 mm	338 mm	10
	Adult	5.0 mm	7.3 mm	16.5 mm	358 mm	10
		5.5 mm	8.0 mm	18.5 mm	368 mm	10
		6.0 mm	8.7 mm	20.5 mm	415 mm	10
		6.5 mm	9.3 mm	24.0 mm	420 mm	10
		7.0 mm	10.0 mm	26.0 mm	430 mm	10
		7.5 mm	10.7 mm	26.0 mm	445 mm	10
		8.0 mm	11.3 mm	28.0 mm	455 mm	10



**PEDI SAFETY SILK**

**RÜSCH**

REF.		ORDER SIZE/I.D.	O.D.	BLACK TIP	LENGTH	QTY
100480	Paediatric	2.0 mm	3.3 mm	20 mm	140 mm	10
		2.5 mm	4.0 mm	20 mm	160 mm	10
		3.0 mm	4.7 mm	20 mm	180 mm	10
		3.5 mm	5.3 mm	20 mm	200 mm	10
		4.0 mm	6.0 mm	30 mm	220 mm	10
		4.5 mm	6.7 mm	30 mm	240 mm	10
		5.0 mm	7.3 mm	40 mm	270 mm	10
		5.5 mm	8.0 mm	40 mm	280 mm	10



**NASOPHARYNGEAL SAFETY SILK**

**RÜSCH**

REF.		ORDER SIZE/I.D.	O.D.	LENGTH	QTY
125410	Paediatric	2.5 mm	3.97 mm	60 mm	10
		3.0 mm	4.64 mm	72 mm	10
		3.5 mm	5.30 mm	85 mm	10
		4.0 mm	5.97 mm	105 mm	10
		4.5 mm	6.29 mm	110 mm	10
	Adult	5.0 mm	6.63 mm	115 mm	10
		5.5 mm	7.29 mm	125 mm	10
		6.0 mm	7.96 mm	130 mm	10
		6.5 mm	8.63 mm	140 mm	10
		7.0 mm	9.29 mm	155 mm	10
		7.5 mm	9.96 mm	165 mm	10
		8.0 mm	10.63 mm	170 mm	10
		8.5 mm	11.29 mm	175 mm	10
		9.0 mm	11.95 mm	180 mm	10



Teleflex is a global provider of medical products designed to enable healthcare providers to protect against infections and improve patient and provider safety. The company specialises in products and services for vascular access, respiratory, general and regional anaesthesia, cardiac care, urology and surgery. Teleflex also provides specialty products for device manufacturers.

The Teleflex family of brands includes ARROW®, BEERE MEDICAL®, DEKNATEL®, GIBECK®, HUDSON RCI®, KMEDIC®, PILLING®, PLEUR-EVAC®, RÜSCH®, SHERIDAN®, SMD™, TAUT®, TFX OEM®, VASONOVA™ and WECK®, all of which are trademarks or registered trademarks of Teleflex Incorporated or its affiliates. Teleflex Incorporated (NYSE:TFX) has annual revenues of approximately \$1.8 billion and serves healthcare providers in 130 countries.

Teleflex global operations: Austria, Belgium, Canada, China, Czech Republic, France, Germany, Greece, India, Ireland, Italy, Japan, Mexico, Netherlands, Portugal, Singapore, Slovak Republic, South Africa, Spain, Switzerland, United Kingdom, Uruguay and USA.

#### **YOUR CONTACTS FOR EUROPE, THE MIDDLE EAST AND AFRICA (EMEA):**

##### **TELEFLEX MEDICAL HEADQUARTERS EMEA, IRELAND**

Teleflex Medical Europe Ltd., IDA Business Park, Athlone, Co Westmeath  
Phone +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73  
orders.intl@teleflexmedical.com

**AUSTRIA** +43 (0)1 402 47 72

**BELGIUM** +32 (0)2 333 24 60

**CZECH REPUBLIC** +420 (0)495 759 111

**FRANCE** +33 (0)5 62 18 79 40

**GERMANY** +49 (0)7151 406 0

**GREECE** +30 210 67 77 717

**ITALY** +39 0362 58 911

**NETHERLANDS** +31 (0)88 00 215 00

**PORTUGAL** +351 22 541 90 85

**SLOVAK REPUBLIC** +421 (0)3377 254 28

**SOUTH AFRICA** +27 (0)11 807 4887

**SPAIN** +34 918 300 451

**SWITZERLAND** +41 (0)31 818 40 90

**UNITED KINGDOM** +44 (0)1494 53 27 61

For detailed information see [www.teleflex.com](http://www.teleflex.com)

The products in this catalogue are only available in the EMEA (Europe, Middle East, Africa) region. For further information contact your local representative. All data current at time of printing (03/2011). Subject to technical changes without further notice.

94 04 56 - 00 00 01 · REV A · MC / WM · 03 11 02