



T-Bag


An oxygen enhancement device designed to increase FiO_2 concentrations post surgery.

- ▶ Single use
- ▶ Lightweight
- ▶ Latex-free
- ▶ Reliable FiO_2 concentrations of 70% at 6 l/min and 50% at 3 l/min
- ▶ Holds approximately 300 ml O_2 in its reservoir
- ▶ No significant rebreathing
- ▶ Confirmation of respiration through T-Bag movement
- ▶ Available in boxes of 100 units






T-BAG - OXYGEN ENHANCEMENT DEVICE

	REF.	QTY
	105230-000010	100

Distributed by:



TELEFLEX HEADQUARTERS INTERNATIONAL, IRELAND
 Teleflex Medical Europe Ltd., IDA Business and Technology Park,
 Dublin Road, Athlone, Co. Westmeath, Ireland
 Phone +353 (0)9 06 46 08 00 Fax +353 (0)14 37 07 73
 orders.intl@teleflex.com



Manufactured by:
 GMMI SDN BHD
 No. 12, Persiaran Industri Rapat 2,
 Perindustrian Ringan Sri Rapat,
 31350 Ipoh, Perak, Malaysia



Consult IFU on this website:
www.LMACO.com

Teleflex www.teleflex.com LMA www.lmaco.com

Copyright © 2015 Teleflex Incorporated. All rights reserved. LMA, LMA Better by Design are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.

94 05 98 – 00 00 01 . REV A · LMA / PDF · 07 15 PDF



Clinical evidence



For the latest clinical evidence on T-Bag Oxygen Enhancement Device

www.lmaco.com/evidence



For the latest digital case reports, educational videos and clinician testimonials

www.youtube.com/LaryngealMaskAirway



For the latest news from LMA

www.facebook.com/LMAInternational



For product information and access to product instructions for use

www.lmaco.com



For information on other products within the Teleflex product portfolio

www.teleflex.com

