



TELEFLEX NON-DEHP PRODUCTS

For aerosol and oxygen therapy

AT A GLANCE

- crush-resistant star-lumen tubing
- identification facilitated by the letter “P” placed at the end of each product reference
- DEHP-free
- 5-year shelf life

WHY PROVIDE NON-DEHP PRODUCTS?

DEHPs are commonly used as plasticizers, particularly in PVC, to provide the desired flexibility, extensibility and elasticity.

In recent years, their frequent use in the production of medical devices has induced responsible international authorities to study their potential toxicity to human health and adapt applicable regulations accordingly:

The revised MDD93/42/EC directive, as amended by the 2007/47/EC directive, includes a new requirement for certain medical devices containing DEHPs.

Our oxygen & aerosol devices meet these requirements, thus enabling you to avoid the risks linked to the use and/or exposure to such substances.

TELEFLEX'S SAFE, NON-DEHP PRODUCTS

To ensure the devices can be used safely in many high-risk situations, Teleflex has developed a range of safe products dedicated to aerosol and oxygen therapy. In addition to reliability and performance levels that are comparable with traditional products, these devices have additional technical characteristics which provide greater safety for the patient and medical staff such as a safe, crush-resistant and 5-year shelf life.

HOW TO IDENTIFY THE PRODUCTS?

Teleflex non-DEHP aerosol and oxygen therapy products are easily identifiable thanks to the letter “P” placed at the end of each product reference*.



View of the inner star-lumen

* please see product codes on the back page of the document

OXYGEN THERAPY

ADULT NASAL CANNULAS

Anatomically-designed, curved, soft, light nasal tips for greater comfort. Star-lumen tubing.

ADULT NASAL CANNULAS

REF.	DESCRIPTION	QTY
1103P	non-flared nasal tip	50
1104P	flared nasal tip	50

MEDIUM AND HIGH CONCENTRATION MASKS

Flexible and transparent oxygen masks with safety vent and 210 cm-long, star-lumen oxygen-supply tubing.

MEDIUM AND HIGH CONCENTRATION MASKS

REF.	DESCRIPTION	QTY
1041P	adult medium concentration mask	50
41042P	paediatric medium concentration mask	50
1059P	adult high concentration mask with check valve and reservoir bag	50
41098P ¹	air entrainment mask with tubing	50

¹Venturi System

OXYGEN SUPPLY TUBING

Star-lumen oxygen supply tubing.

OXYGEN SUPPLY TUBING

REF.	DESCRIPTION	QTY
41113P	210 cm	50
1115P ²	210 cm	50

²Plain oxygen supply tubing

AEROSOL THERAPY

NEB EASY SMALL VOLUME NEBULISER

Small volume nebuliser (6 ml reservoir) 200 cm-long star-lumen tubing.

NEB EASY SMALL VOLUME NEBULISER

REF.	DESCRIPTION	QTY
41892P	with T-shaped mouthpiece	50
41893P	with adult mask	50
41894P	with paediatric mask	50

ADULT AEROSOL MASK

Anatomically designed flexible mask with chin cup and nose clip. Connector accepts 22 mm aerosol tubing.

ADULT AEROSOL MASK

REF.	DESCRIPTION	QTY
1083P	adult mask	50

For more information about our products, please contact our customer service department.

DISTRIBUTED BY:

TELEFLEX MEDICAL HEADQUARTER INTERNATIONAL, IRELAND

IDA Business & Technology Park · Dublin Road · Athlone · Co Westmeath
Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com

UNITED KINGDOM Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com

SOUTH AFRICA Tel. +27 (0)11 807 4887 · assist.africa@teleflex.com

WWW.TELEFLEX.COM

94 06 67 · 00 00 01 · REV A · 01 13 01 · All data current at time of printing. Subject to technical changes without further notice.



Ref. 1059P

Ref. 41113P



Ref. 41892P



Ref. 1083P