



Arrow  
**Midline**  
 Simple and Safe Peripheral Venous Access

It is the ambition of Teleflex to always be one step ahead. Our products are designed, manufactured and distributed for healthcare customers. Striving for precision in detail to provide safety coupled with the comfort for both patient and caregiver is our mission.

**Stay ahead of risks**

Inspired by the same concern with safety that has made our central venous catheter kits the industry leader, the Arrow Midline kit provides you with a comprehensive collection of high-quality components.

**Simplicity of use for maximum quality**

With this in mind, Teleflex is glad to present the newly developed Basic Kit including the Arrow Midline. The kit only contains components of high quality.

The new Arrow GlideThru peelable sheath, featuring the new twist lock connection and the new hub design, helps assist in a smooth transition from dilator to sheath, thus reducing the need for a skin nick.

The sharpsaway II locking disposal protects the operator from needlestick injuries.

These are only an example of the attention that Teleflex puts in to designing its kits to make the lives of clinicians easier.



Benefits




**The Clinician**  
 the IV therapy solution can be tolerated by a larger peripheral vessel



**Your Institution**  
 therapy will last no longer than 1-4 weeks



**The Patient**  
 helps to protect from multiple patient needle stick



## Kit contents

- |   |  |  |
|---|--|--|
| <ul style="list-style-type: none"> <li>① Nitinol spring wire assembly .018" x 45cm<br/>atraumatic and kink resistant</li> <li>② GlideThru peelable sheath<br/>smooth dilator-to-sheath transition</li> <li>③ Echogenic insertion needle, 7cm<br/>enhanced visibility for ultrasound guided<br/>insertion</li> </ul> | <ul style="list-style-type: none"> <li>④ Needle-free Connectors<br/>reduce the risk of contamination and blood<br/>reflux</li> <li>⑤ SharpsAway II locking disposal<br/>protects the operator from needlestick injuries</li> <li>⑥ 10ml syringe</li> <li>⑦ Safety scalpel #11</li> </ul> | <ul style="list-style-type: none"> <li>⑧ Arrow Midline catheter</li> <li>⑨ Sutureless fixation device<br/>provides atraumatic and stable fixation</li> <li>⑩ Guillotine catheter trimmer</li> <li>⑪ Tourniquet and tape measure</li> </ul> |
|---|--|--|

## Ordering information

ARROW

REF.	N° OF LUMENS	OUTER DIAMETER	INNER DIAMETER	FLOW RATE*	LENGTH	QTY
EU-02031-ML	1	3 Fr.	20 G	300 mL/hr	20 cm	5
EU-02041-ML	1	4 Fr.	18 G	1690 mL/hr	20 cm	5
EU-02051-ML	1	5 Fr.	18 G	3060 mL/hr	20 cm	5

\* Gravity Flow Rate



teleflex.com

### References:

- Grip-Lok is a registered trademark of Zefon International, Inc.
- Trerotola S., Stravropoulos S., Mondeschein J., Patel A., Fishman N., Fuchs B., Kolansky D., Kasner S., Pryor J., Chittams J.
- "Triple-Lumen Peripherally Inserted Central Catheter in Patients in the Critical Care Unit: Prospective Evaluation." Radiology, 2010; 256 (1.2).
- MicroCLAVE Neutral Displacement Connectors is an exclusive proprietary product manufactured by ICU Medical

Teleflex, the Teleflex logo, and Arrow are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this catalogue may not be available in all countries. Please contact your local representative. All data current at time of printing (10/2016). Subject to technical changes without further notice. © 2016 Teleflex Incorporated. All rights reserved. 94 11 36 - 00 00 01 - REV B · 10 16 PDF

**Teleflex**

Distributed by:  
**Teleflex Headquarters International, Ireland** · Teleflex Medical Europe Ltd. · IDA Business & Technology Park  
 Dublin Road · Athlone · Co Westmeath · Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com  
**United Kingdom** Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com  
**South Africa** Tel. +27 (0)11 807 4887 · assist.africa@teleflex.com