



ARROW MIDLINE

A customised solution to meet clinicians' needs

AT A GLANCE

THE ARROW MIDLINE

- made of medical grade, flexible polyurethane
- clearly marked as a midline on the extension line
- now available in MSB (maximal sterile barrier) configuration
- the safety hypo needles helps to protect from multiple patient needle stick

KEY BENEFITS

Midlines are appropriate when

- facing difficult IV catheter insertions
- therapy will last no longer than 1-4 weeks
- the IV therapy solution can be tolerated by a larger peripheral vessel

COMPREHENSIVE AND EASY TO USE FOR THE INSERTER

It is the ambition of Teleflex to always be one step ahead. Our products are designed, manufactured and distributed for healthcare customers. Striving for precision in detail to provide safety coupled with the comfort for both patient and caregiver is our mission.

The Arrow Midline kit includes virtually everything you need to perform a midline placement – while minimising the risk to your patient and yourself.

STAY AHEAD OF RISKS

Inspired by the same concern with safety that has made our central venous catheter kits the industry leader, the Arrow Midline MSB kit provides you with a comprehensive collection of high-quality components designed to help minimize the risk of catheter-related infection and sharps-related injuries. It also puts nearly everything you need for Midline insertion and securement right at your fingertips.

PROVEN PROTECTION THANKS TO MAXIMAL BARRIER PRECAUTIONS

According to the Institute for Healthcare Improvement, "... two studies of central line insertions have shown that the application of maximal barrier precautions substantially reduces the odds of developing a bloodstream infection. Mermel et al. demonstrated that the odds ratio was 2,2 times greater for infection without maximal barrier precautions, while Raad et al. demonstrated a 6,3 times greater likelihood for infection without precautions."¹





KIT CONTENTS:

- ① Nitinol spring wire assembly .018" (45 cm) atraumatic and kink resistant
- ② GlideThru® peelable sheath smooth dilator-to-sheath transition
- ③ Echogenic insertion needle, 7cm enhanced visibility for ultrasound guided insertion
- ④ Insertion needle, 4cm
- ⑤ Needle-free Connectors reduce the risk of contamination and blood reflux
- ⑥ SharpsAway II™ locking disposal protects the operator from needlestick injuries
- ⑦ Safety injection needles protecting the operator from needlestick injuries
- ⑧ Filter needle
- ⑨ 2ml syringe
- ⑩ 5ml syringe
- ⑪ 10ml syringe
- ⑫ Safety scalpel #11
- ⑬ Secondsite™ adjustable hub fastens catheter at insertion site anywhere along the catheter body
- ⑭ Arrow Midline catheter 4Fr x 20cm (Polyurethane)
- ⑮ Transparent dressing – 10cm x 15.5cm
- ⑯ Sutureless fixation device provides atraumatic and stable fixation
- ⑰ Dressing – 6.5cm x 5cm
- ⑱ Maximal barrier tools: mask with eye shield, surgical gown, exclusive easy-to-use head-to-toe drape, fenestrated drape
- ⑲ Sterile ultrasound probe cover + gel
- ⑳ Sterile gauzes

ARROW MIDLINE MSB KIT

ARROW

REF.	DESCRIPTION	FRENCH (FR)	LENGTH (CM)	QTY
UK-00420-MID	Arrow Midline maximal sterile barrier kit	4	20	4

References:

¹ "Getting Started Kit: Prevent Central Line Infections." The Institute for Healthcare Improvement, 2005, p. 5.

- Grip-Lok® is a registered trademark of Zefon International, Inc.
- SafetyEclipse™ is a registered trademark of Becton, Dickinson and Company.
- Tegaderm® is a registered trademark of 3M Company.
- MicroCLAVE® Neutral Displacement Connectors is an exclusive proprietary product manufactured by ICU Medical

Teleflex and Arrow are trademarks of Teleflex Incorporated or its affiliates registered in the U.S. and/or other countries.
© 2015 Teleflex Incorporated. All rights reserved

DISTRIBUTED BY:

TELEFLEX HEADQUARTERS INTERNATIONAL, IRELAND · Teleflex Medical Europe Ltd. · IDA Business & Technology Park
Dublin Road · Athlone · Co Westmeath · Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com

UNITED KINGDOM Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com

SOUTH AFRICA Tel. +27 (0)11 807 4887 · assist.africa@teleflex.com

WWW.TELEFLEX.COM

94 11 36 - 00 UK 01 · REV A · 04 15 0.2 · All data current at time of printing. Subject to technical changes without further notice.

Teleflex