



ARROW® OnControl®
POWERED BONE ACCESS



BONE MARROW BIOPSIES

The new standard starts here

Teleflex

ARE INCONSISTENT SPECIMEN SIZE AND MULTIPLE PROCEDURES COMMON IN YOUR FACILITY?

Trusted brands make Teleflex a reliable and strong partner. Built on a solid tradition of innovation, Teleflex is a global leader in superior medical supplies designed to help providers minimise risk and maximise outcomes for their patients. Our understanding of the importance to our customers of a full range of products has led to the development of a unique line of medical devices, all of which complement one another.


WE'RE WITH YOU EVERY STEP OF THE WAY.

The ARROW® OnControl® powered bone marrow biopsy system, supported by our experienced clinical training staff, has successfully raised the standard for bone marrow biopsies and aspirations in leading cancer centers, while greatly improving the patient experience.¹



You can find further details and the technical specifications of our products within this brochure.

TELEFLEX – HIGH QUALITY MEDICAL SUPPLIES FROM A SINGLE SOURCE

CORE SPECIMEN SIZE BY VOLUME (MM³)



BONE MARROW BIOPSY SPECIMENS



Bone marrow biopsy specimens were obtained using the ARROW® OnControl® powered bone marrow biopsy system 11g needle. The specimens measured 1.7 cm in length.

Individual results may vary


- more consistent, high quality specimens
- longer, wider, and larger volume specimens^{1,2,3,4*}
- more usable area for diagnosis^{1*}

NOW YOU CAN OBTAIN CONSISTENTLY LARGER, HIGH QUALITY CORE SPECIMENS

ASSISTS CLINICIANS IN MAKING THE DIAGNOSIS

Unlike specimens obtained from manual biopsy needles, specimens provided by the ARROW® OnControl® Powered Bone Marrow Biopsy System are high quality and consistently larger in size,^{1,2,3,4} providing more usable area for diagnosis.¹ These longer, wider, and greater volume specimens^{1,2,3,4} reduce second attempt procedures required^{1,4} when the specimen fails to capture.

Clinicians can be confident in treatment plans based on a diagnosis specimen for patients with blood-borne cancers.



ARROW® OnControl® Bone Marrow Needle:
specially designed threaded cannula captures
and holds core specimens

ENSURES BETTER OVERALL CLINICIAN AND PATIENT EXPERIENCE

Patients with blood-borne cancer frequently describe the bone marrow biopsy and aspiration as the “worst part of their treatment.” The ARROW® OnControl® Powered Bone Marrow Biopsy System helps improve patient satisfaction⁴ and promote compliance with ongoing and future testing, by significantly reducing the pain during and after the procedure,^{3,4} and reducing the need for “second attempts” that often occur when a specimen is of insufficient size or fails to capture^{1,3}.


ARROW® OnControl® Powered Bone Marrow Biopsy System has a faster procedure time and increased control, as compared to manual biopsy needles, reduces the physical requirements and effort needed by clinicians performing the biopsy, even with hard bone.

- patented threaded cannula design grabs and holds core specimens from soft bone
- increased user control and reduced physical requirements to obtain specimens, even with hard bone
- reduced insertion pain,¹ less overall pain,³ and helps promote patient compliance^{3,4}
- greater overall patient satisfaction^{4*}

IMPROVES CAREGIVER EFFICIENCY

When compared to biopsies performed with manual needles, the ARROW® OnControl® Powered Bone Marrow Biopsy System provides improved specimen quality and shortened procedure time, that may help clinicians perform a greater number of biopsies per day for patients with blood-borne cancers and benign tumors.

ARROW® OnControl®
Core Specimen**



Manual Needle
Core Specimen


* Compared to manual bone marrow biopsy procedures

* Both specimen sizes from same healthy subject, same provider. Specimen sizes are most representative of specimen means in a published healthy subject study.


- fewer second-attempt procedures required^{3*}
- 55% faster procedure time to improve efficiency^{1*}
- easy for clinicians to use, regardless of physical strength



A complete solution for powered Bone Marrow Biopsies

ONCONTROL POWER DRIVER




ONCONTROL POWER DRIVER		ARROW
	REF.	QTY
	9401	1

BONE MARROW BIOPSY SYSTEM

BIOPSY SYSTEM TRAYS, 11 G			ARROW
	REF.	LENGTH	QTY
	9408-EU-006	102 mm	6
	9411-EU-006	152 mm	6

BIOPSY SYSTEM NEEDLE SETS, 11 G			ARROW
	REF.	LENGTH	QTY
	9407-EU-006	102 mm	6
	9402-EU-006	152 mm	6

BONE MARROW ASPIRATION SYSTEM

ASPIRATION SYSTEM NEEDLE SETS, 15 G			ARROW
	REF.	LENGTH	QTY
	9425-VC-006	25 mm	6
	9468-VC-006	68 mm	6
	9490-VC-006	90 mm	6

Some products not available in all markets, subject to availability. Contact your local representative.

ARROW® OnControl® powered bone access systems should only be used by clinicians familiar with the complications, limitations, indications, and contraindications of the indicated procedures.

- Swords RT, Anguita J, Higgins RA, et al.** A prospective randomized study of a rotary powered device (OnControl) for bone marrow aspiration and biopsy. *J Clin Pathol* 2011;64(9):809-13. doi:10.1136/jclinpath-2011-200047.*
- Berenson JR, Yellin O, Blumenstein B, et al.** Using a powered bone marrow biopsy system results in shorter procedures, causes less residual pain to adult patients, and yields larger specimens. *Diagnostic Pathology* 2011;6:23.*
- Miller LJ, Philbeck TE, Montez DF, et al.** Powered bone marrow biopsy procedures produce larger core specimens, with less pain, in less time than with standard manual devices. *Hematology Reports* 2011;3:e8.t.*
- Reed LJ, Raghupathy R, Strakhan M, et al.** The OnControl bone marrow biopsy technique is superior to the standard manual technique for hematologists-in-training: a prospective, randomized comparison. *Hematology Reports* 2011;3(e21). doi:10.4081/hr.2011.e21.*

* Research sponsored by Teleflex Incorporated or its affiliates, including Vidacare LLC.

For more information, visit us at: arrowoncontrol.com.

Teleflex, ARROW and OnControl are trademarks or registered trademarks of Teleflex Incorporated or its affiliates registered in the U.S. and/or other countries.

TELEFLEX HEADQUARTERS INTERNATIONAL, IRELAND

Teleflex Medical Europe Ltd., IDA Business and Technology Park,
Dublin Road, Athlone, Co Westmeath

Phone +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73

orders.intl@teleflex.com

AUSTRALIA/NEW ZEALAND +61 (0)3 9081 0600

AUSTRIA +43 (0)1 402 47 72

BELGIUM +32 (0)2 333 24 60

CHINA (SHANGHAI) +86 (0)21 6163 0965

CHINA (BEIJING) +86 (0)10 6418 5699

CZECH REPUBLIC +420 (0)495 759 111

FRANCE +33 (0)5 62 18 79 40

GERMANY +49 (0)7151 406 0

GREECE +30 210 67 77 717

INDIA +91 (0)44-2836 5040

ITALY +39 0362 58 911

JAPAN +81 (0)3 3379 1511

KOREA +82 2 536 7550

NETHERLANDS +31 (0)88 00 215 00

PORTUGAL +351 22 541 90 85

SINGAPORE +65 6439 3000

SLOVAK REPUBLIC +421 (0)3377 254 28

SOUTH AFRICA +27 (0)11 807 4887

SPAIN +34 918 300 451

SWITZERLAND +41 (0)31 818 40 90

UNITED KINGDOM +44 (0)1494 53 27 61

For detailed information see www.arrowoncontrol.com

The products in this catalogue may not be available in all countries.
Please contact your local representative. All data current at time of
printing (09/2015). Subject to technical changes without further notice.

940785 - 00 00 01 · REV A · MC / SF · 07 15 03

Teleflex