



**ARROW® EZ-IO®**  
INTRAOSSUEOUS VASCULAR ACCESS



## EZ-IO® NEEDLE + STABILIZER KIT

### AT A GLANCE

Easy-to-use kit includes:

- 15 G sterile EZ-IO needle set (45 mm, 25 mm, and 15 mm)
- EZ-Stabilizer dressing
- EZ-Connect® extension set
- NeedleWISE® 1-port
- patient-wristband
- same packaging size as the current EZ-Stabilizer dressing



M-872 Rev A

### STREAMLINES AND SIMPLIFIES THE EZ-IO VASCULAR APPLICATION PROCESS

The ARROW EZ-IO Needle + Stabilizer Kit provided by Teleflex, is a single easy-to-use kit that streamlines and simplifies the EZ-IO vascular access application process and eliminates the manual securement of the EZ-IO needle.

Each kit comes with an EZ-IO needle set and an EZ-Stabilizer® dressing. The EZ-Stabilizer dressing provides secure placement of the EZ-IO intraosseous vascular access needle, preventing potential dislodgement due to accidental impact on the needle or other external factors. It can be placed regardless of location and is recommended for use with all EZ-IO intraosseous vascular access needle placements.

### KEY BENEFITS

- provides secure placement of the EZ-IO needle
- eliminates the need for manual securement of the IO needle
- reduces the number of products for clinicians to carry for IO procedures
- streamlines the EZ-IO vascular access system application process
- saves clinician's time and maintains consistency of application

EZ-IO NEEDLE + STABILIZER KIT			ARROW
	REF.	DESCRIPTION	QTY
	9079P	45 mm needle + stabilizer	5
	9001P	25 mm needle + stabilizer	5
	9018P	15 mm needle + stabilizer	5

For more information, please visit [www.arrowezio.com](http://www.arrowezio.com). Teleflex and EZ-IO are trademarks or registered trademarks of Teleflex Incorporated or its affiliates. © 2014 Teleflex Incorporated.

**DISTRIBUTED BY:**  
**TELEFLEX HEADQUARTERS INTERNATIONAL, IRELAND** · Teleflex Medical Europe Ltd. · IDA Business & Technology Park  
 Dublin Road · Athlone · Co Westmeath · Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · [orders.intl@teleflex.com](mailto:orders.intl@teleflex.com)  
**UNITED KINGDOM** Tel. +44 (0)14 94 53 27 61 · [info.uk@teleflex.com](mailto:info.uk@teleflex.com)  
**SOUTH AFRICA** Tel. +27 (0)11 807 4887 · [assist.africa@teleflex.com](mailto:assist.africa@teleflex.com)  
**WWW.TELEFLEX.COM**  
 94 07 75 - 00 00 01 · REV A · 09 14 05 · All data current at time of printing. Subject to technical changes without further notice.

