



GIBECK ISO-GARD® FILTER

Hydrophobic bacterial/viral filters



AT A GLANCE

Iso-Gard filter,
Iso-Gard filter small

- filtration efficiency (99.9999%)
- for tidal volumes of 150–1000 ml
- with CO₂ port
- small, compact, easy to handle
- available in straight and angled version
- choice for adults and children
- latex-free

Iso-Gard Filters are also available in several options with different Flex tubes.

PROTECT THE PATIENT, MEDICAL PERSONNEL AND VENTILATION EQUIPMENT FROM CONTAMINATION.

The high filtration efficiency (99.9999%) of Teleflex Gibeck Iso-Gard filters has been certified by independent institutes who have carried out extensive tests.

The filters, available in two versions, named Iso-Gard filter and Iso-Gard filter small, comply with all the required quality standards. They are ideal for patient- and machine close application in anaesthesia and for machine close application in the ICU.

The Iso-Gard filters are extremely lightweight and easy to handle.

Thanks to the low dead space the Iso-Gard filter small is suitable for both children and adults.



PRODUCT INFORMATION

Gibeck

	ISO-GARD® FILTER		ISO-GARD® FILTER SMALL	
	ANGLED	STRAIGHT	ANGLED	STRAIGHT
Ref. sterile clean	18211 18212	19211 19212	18511 18512	19511 19512
Dead space	30ml	26ml	21ml	20ml
Weight	23g	22g	16g	15g
Connectors	22 M/15 F-15 M/22 F		22 M/15 F-15 M	
Medium	hydrophobic depth filter made of polypropylene, electrostatically charged		hydrophobic depth filter made of polypropylene, electrostatically charged	
Tidal volume	150–1000 ml		150–1000 ml	
Resistance to flow	at 60l/min., 1.6 cm H ₂ O		at 30l/min., 1.9 cm H ₂ O	
Bacterial efficiency	> 99.9999 %		> 99.9999 %	
Viral efficiency	> 99.999 %		> 99.999 %	
Qty (box/case)	250/25		200/50	

DISTRIBUTED BY:

TELEFLEX MEDICAL HEADQUARTER EMEA, IRELAND · IDA Business & Technology Park · Dublin Road · Athlone · Co Westmeath

Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com

UNITED KINGDOM Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com

SOUTH AFRICA Tel. +27 (0)11 807 4887 · assist.africa@teleflex.com

WWW.TELEFLEX.COM

94 03 72 - 00 00 01 · REV A · 01 13 01· All data current at time of printing. Subject to technical changes without further notice.

teleflex