



Arrow
Sharps Away II
 Locking Disposal Cup

Daily Risk of a Healthcare Worker

Everyday work of healthcare staff puts them at risk of serious infections, with more than 30 potentially dangerous pathogens, including hepatitis B, hepatitis C and HIV, as a result of needlestick injuries¹. Needlestick injuries are a very serious occupational hazard for healthcare workers.

More than one million needlestick injuries are estimated to occur in the European Union each year. Needlestick injuries are one of the most common and serious risks to healthcare workers in Europe and represent a high cost for health systems and society in general.²

Easy and Safe Protection

Arrow helps prevent infection and reduce risk of needle stick injuries with the Arrow Sharps Away II Locking Disposal Cup which complies with the latest EU Directive 2010/32/EU-prevention from sharp injuries in the hospital and healthcare sector.

The Arrow Sharps Away II Locking Disposal Cup provides an easy and safe way of moving used sharps from procedure to sharps disposal container. It allows safe disposal of used sharps to minimise the risk of exposure.

Benefits



The Clinician
 Prevent risk of needle stick injuries



Your Institution
 Comply with the latest standards

Product Details

The Sharps Away II Disposal Cup has an adhesive surface on the bottom side which allows a fixed position for a one-handed technique.

Once placed into disposal cup, needles will be automatically secured in place so that they cannot be reused.

- Each section can be used for one needle or canula (maximum 8 pcs. per one Sharps Away II).
- 15 G - 30 G



Easily move used sharps to disposal container while reducing exposure to sharps injuries

Sharps Away II

REF.	DESCRIPTION	QTY
SA-00011	Sharps Away II Locking Disposal Cup	250 pcs.

¹ http://www.europeanbiosafetynetwork.eu/EBN_Implementation_Update.pdf

² Prevention of Sharps Injuries in the Hospital and Healthcare Sector/ Implementation Guidance for the EU Framework Agreement, Council Directive and Associated National Legislation



teleflex.com

Teleflex, the Teleflex logo, Arrow and Sharps Away are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this catalogue may not be available in all countries. Please contact your local representative. All data current at time of printing (03/2017). Subject to technical changes without further notice.

© 2017 Teleflex Incorporated. All rights reserved. 94 11 38 - 00 00 01 - REV A - 03 17 PDF

Teleflex

Distributed by:

Teleflex Headquarters International, Ireland · Teleflex Medical Europe Ltd. · IDA Business & Technology Park

Dublin Road · Athlone · Co Westmeath · Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com

United Kingdom Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com

South Africa Tel. +27 (0)11 807 4887 · assist.africa@teleflex.com