



LMA Classic™

The classic laryngeal mask.

LMA Classic™ was the original laryngeal mask released into the market. LMA Classic™ is a first generation, versatile, re-usable airway with aperture bars used in millions of surgical cases each year.

- ▶ Eight sizes, from neonate to large adult
- ▶ Not made with natural rubber latex and re-usable up to 40 times
- ▶ Aperture bars designed to prevent the blockage of the airway by the epiglottis¹
- ▶ Soft, silicone cuff
- ▶ Facilitates smooth emergence from anaesthesia²
- ▶ Minimal haemodynamic response²
- ▶ Part of the ASA difficult airway algorithm³


Consider using LMA Classic™ for:

- ▶ Theatre
- ▶ Crash carts
- ▶ Emergency departments
- ▶ Pre-hospital
- ▶ Ambulance



“Advantages over the TT included: increased speed and ease of placement by inexperienced personnel; increased speed of placement by anaesthetists; improved haemodynamic stability at induction and during emergence; minimal increase in intraocular pressure following insertion; reduced anaesthetic requirements for airway tolerance; lower frequency of coughing during emergence; improved oxygen saturation during emergence; and lower incidence of sore throat in adults.”

Brimacombe J., 1995




LMA Classic™: Product specification

Mask size	Product code	Patient size	Maximum cuff volume (air)*
1	100010	Neonates/infants up to 5 kg	4 ml
1.5	100015	Infants 5-10 kg	7 ml
2	100020	Infants/children 10-20 kg	10 ml
2.5	100025	Children 20-30 kg	14 ml
3	100030	Children 30-50 kg	20 ml
4	100040	Adults 50-70 kg	30 ml
5	100050	Adults 70-100 kg	40 ml
6	100060	Large adults over 100 kg	50 ml

*These are maximum volumes that should never be exceeded. It is recommended that the cuff be inflated to a maximum of 60 cm H₂O intracuff pressure.

Distributed by:



PAA-2500-005 Rev B IW 201402

TELEFLEX HEADQUARTERS INTERNATIONAL, IRELAND
 Teleflex Medical Europe Ltd., IDA Business and Technology Park,
 Dublin Road, Athlone, Co Westmeath, Ireland
 Phone +353 (0)9 06 46 08 00 Fax +353 (0)14 37 07 73
 orders.intl@teleflex.com



Manufactured by:
 The Laryngeal Mask Company Limited
 Le Rocher, Victoria, Mahé, Seychelles



Consult IFU on this website:
www.LMACO.com



Printed on paper from sustainable sources

Teleflex www.teleflex.com LMA www.lmaco.com

Copyright © 2014 Teleflex Incorporated. All rights reserved. LMA, LMA Classic, LMA Better by Design are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.

94 06 91 - 00 00 01 - REV B - LMA / PDF - 03 14 PDF

References:

1. Benumof J.L. CJA 2003; **50**: 968.
2. Brimacombe J. CJA 1995; **42**: 1017-1023.
3. American Society of Anaesthesiologists Task Force on Management of the Difficult Airway. Anesthesiol 2003; **98**: 1269-1277.



Clinical evidence



For the latest clinical evidence on LMA Classic™

www.lmaco.com/evidence



For the latest educational videos and clinician testimonials

www.youtube.com/LaryngealMaskAirway



For the latest news from LMA

www.facebook.com/LMAInternational



Better by Design



For product information and access to product instructions for use

www.lmaco.com



For information on other products within the Teleflex product portfolio

www.teleflex.com

