



Rüsch
PercuTwist
One-Step Dilator

Safe, Intuitive and Cost Effective

The Rüsch PercuTwist One-Step Dilator from Teleflex is a patented single-step dilator that has made percutaneous tracheostomy a controllable and genuinely safe procedure and has thus set new standards.

It comes with an application which is surprisingly simple and gentle. Enhanced cosmetic results are often achieved in comparison with a surgical procedure¹, as the access site requires an incision of less than 2 cm. Its innovative design with a dilating thread makes it possible to lift the tracheal anterior wall while simultaneously dilating through simple rotation. The hydrophilic coating reduces friction and allows smoother access into the trachea.

Available in a choice of 3 sizes to enable the device to be used most effectively for a wide range of adult patients.

An unrestricted endoscopic view of the dilation site is guaranteed at any given time.

Free Choice

The PercuTwist SOLO Set allows you choose the tube to meet your patient's clinical need. The dilators have been designed to mirror the O.D. of the chosen tube, to facilitate smooth and gentle placement in the stoma.

Benefits



The Clinician

Intuitive to use & available for a wide variety of tracheostomy tube types.



Your Institution

Can be performed bedside in acute care - reduces the costs associated with moving patients to the OR (for indicated patients, see IFU).



The Patient

Minimal scarring for patients not requiring a long term stoma.

PercuTwist One-Step Dilator Fast Facts

Innovative design

Hydrophilic coating

Three different diameters

Available with Rüscher PercuQuick or CrystalClear Tracheostomy Tube and as an independent kit

Complete PercuTwist Set with PercuQuick Tracheostomy Tube

REF.	SIZE IN MM	QTY
121555	7.0/8.0/9.0	1

Set Components: PercuTwist Set with PercuQuick Tracheostomy Tube

- PercuTwist Dilator
- PercuQuick Tracheostomy Tube with introduction dilator
- Guide wire with J-tip and insertion aid
- Puncturing cannula 17 G x 70 mm/14 G x 52 mm
- Scalpel
- Luer-lock syringe
- Neck band
- Latex-free, sterile, single use



REF. 121555

Complete PercuTwist Set with CrystalClear Tracheostomy Tube

REF.	SIZE IN MM	QTY
121556	7.0/8.0/9.0	1

Set Components: PercuTwist Set with CrystalClear Tracheostomy Tube

- PercuTwist Dilator
- CrystalClear Tracheostomy Tube with introduction dilator
- Guide wire with J-tip and insertion aid
- Puncturing cannula 17 G x 70 mm/14 G x 52mm
- Scalpel
- Luer-lock syringe
- Neck band
- Latex-free, sterile, single use



REF. 121556

Combination Possibilities PercuTwist SOLO Set

REF.	SIZE IN MM	QTY
121557	7.0/8.0/9.0	1
121550	7.0/8.0/9.0	1
120348	7.0/8.0/9.0	1
121501	7.0/8.0/9.0	1
121502	7.0/8.0/9.0	1

Set Components: PercuTwist SOLO Set

- PercuTwist Dilator
- Guide wire with J-tip and insertion aid
- Puncturing cannula
- Scalpel
- Luer-lock syringe
- Neck band
- Latex-free, sterile, single use

The sizes indicated always correspond to the size of the tracheostomy tube used.

References:

1. Silvester et al. Percutaneous versus surgical tracheostomy: A randomized controlled study with long-term follow-up. Crit Care Med. 2006; Aug; 34(8):2145-52.

Teleflex, the Teleflex logo, CrystalClear, Rüscher, PercuQuick and PercuTwist are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this catalogue may not be available in all countries. Please contact your local representative. All data current at time of printing (07/2016). Subject to technical changes without further notice.

© 2016 Teleflex Incorporated. All rights reserved. 94 30 13 - AP 00 01 - REV C - 07 16 PDF

Teleflex

Distributed by:
 Teleflex Headquarters International, Ireland · Teleflex Medical Europe Ltd. · IDA Business & Technology Park
 Dublin Road · Athlone · Co Westmeath · Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com
 Singapore (SEA non-direct sales countries) · Tel. +65 6439 3000 · singaporecustomerservice@teleflex.com
 www.teleflex.com



REF. 121550

REF. 120348

REF. 121501

REF. 121502



REF. 121557