



VPS Rhythm™ DLX Device

PICC placement for those who **don't like surprises**

The VPS Rhythm™ DLX Device is available in 5 configurations

- ✓ Ultrasound with TipTracker™ Technology and ECG
- ✓ Ultrasound with ECG
- ✓ TipTracker™ Technology with ECG
- ✓ ECG only
- ✓ Ultrasound only



The VPS Rhythm™ DLX Device delivers precision tracking and locating with a full breadth of features that help enhance user experience, including:

- **Vascular access ultrasound** that is easy to use, intuitive to navigate, and includes a catheter-to-vessel ratio (CVR) tool²
- **Electromagnetic navigation** that facilitates visualization of the PICC tip as it advances through the vasculature
- **Electrocardiogram (ECG) tip location technology** indicated as an alternative to a confirmatory chest x-ray¹
- **Versatile, portable, and easy to use; designed for bedside placement** of vascular access catheters²

References:

1. RHY-200-OPM-US: VPS Rhythm™ DLX Operator's Manual.
2. Formative evaluation for ultrasound imaging comparison and usability. D023327. July 11th, 2018.

Federal Law (USA) restricts these devices to sale by or on the order of a physician.

Not all products may be available in all countries. For product information please contact la.cs@teleflex.com or your local representative. Refer to the applicable Instructions for Use for the indications approved in your geography. Please check your local regulatory approval status.

Teleflex, the Teleflex logo, Arrow and VPS Rhythm are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 12/2024. © 2024 Teleflex Incorporated. All rights reserved. MC-010355 LA EN

