

HUDSON RCI® GAS SAMPLING LINES

SAMPLING LINES

- PVC sampling lines for measurement of end-tidal CO₂
- Choose between PVC or co-extruded PVC/polyethylene sampling lines for precision measurement of anesthetic gases
- All gas sampling lines include a male luer connector on one end for connection to the patient circuit gas sampling port
- Choice of male or female luer for machine fitting

ITEM NUMBER	LENGTH	INTERNAL DIAMETER	MALE/MALE	FEMALE/MALE	PVC	CO-EXTRUDED
1232H	120"	.050		●	●	
1233H	120"	.050	●		●	
1234H	120"	.047		●		●
1235	120"	.047	●			●
126263	96"	.050	●		●	
126463	96"	.047	●			●
126613	120"	.040	●		●	
126664	96"	.040		●	●	
126713	120"	.060	●		●	