





MICRO ELITE™ SNARE

**The ONLY 0.014" compatible  
guidewire-based snare**

## SECURE CAPTURE. 014" COMPATIBILITY.

The Micro Elite snare features an ultra-small 0.014" profile that permits delivery through microcatheters and balloon catheters. The preassembled, one-piece construction, provides durability while allowing rapid deployment.




# MICRO Elite™ snare

The MICRO Elite Snare is intended for use in the retrieval and manipulation of atraumatic foreign bodies located in the coronary and peripheral cardiovascular system and the extra-cranial neurovascular anatomy.

Model	Loop Diameter	Working Length	O.D.
8300	2mm	180cm	0.014"
8301	4mm	180cm	0.014"
8302	7mm	180cm	0.014"

\*Packaged in quantities of 1 unit per box.



Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Elite is a trademark of Vascular Solutions, Inc.

©2015 Vascular Solutions, Inc. All rights reserved. ML1901 Rev. E 07/15

**This literature is intended for distribution only in the U.S.**



**Vascular Solutions, Inc.**  
6464 Sycamore Court North  
Minneapolis, Minnesota 55369 USA

Customer Service:  
United States: 888.240.6001  
customerservice@vasc.com  
www.vasc.com