


Insighters® Video Laryngoscope

Envision Versatility

Difficult tracheal intubation during direct laryngoscopy has been shown to occur in 1.5–8.5%¹ of patients, while intubation failure has been shown to occur in approximately 1:2000 patients.² Studies show the use of video laryngoscopy devices improves

the glottic view and may reduce airway trauma and the number of failed intubations.³ The Insighters® Video Laryngoscope facilitates intubation while allowing clinicians to share and record real-time high-resolution images and video.



Benefits



The Clinician
Provides versatility to support clinician intubation needs



Your Institution
Video laryngoscopy may lower total costs, reduce length of stay, and decrease the likelihood of postoperative ICU admission⁴



The Patient
Facilitates intubation success³

Insighters® Video Laryngoscope

Technical Specifications

COMPONENT	SPECIFICATIONS
VL Handle	Camera Resolution: 1,000,000 pixels (1280 x 720) Field of View: ≥60° Illumination: LED light source
Sheath	Material: Polycarbonate Sterile, single-use
Small Screen	Display: 3.5-inch LCD medical resistance touch screen, HD color Resolution: 640 x 480 Rotating Angle: Forward 60° / Backward 90° / Leftward 90° / Rightward 180° Built in Battery: Rechargeable lithium battery, working time 4 hours fully charged Battery Recharging Time: Up to 4 hours
Large Screen	Display: 13.3-inch, TFT-LCD color Resolution: 1920 x 1080 Connection Port: Type-C port, HDMI port, headphone port, power socket Built in Battery: Rechargeable lithium battery, working time 6 hours fully charged Battery Recharging Time: Up to 6.5 hours
Cable	Length: 2 meters
NFC Module	Connection Type: Near Field Communication (NFC) wireless Connection Distance: ≤20 meters Built in Battery: Rechargeable lithium battery, working time 4 hours fully charged Battery Recharging Time: Up to 4 hours

Ordering Information

ITEM NUMBER	ITEM DESCRIPTION	PATIENT WEIGHT	ORDER QTY
ARJ12001RU	Handheld Video Laryngoscope		1/EA
ARJ0100	Video Laryngoscope Sheath (SS)	0.7 – 8 kg	20/BX
ARJ0110	Video Laryngoscope Sheath (S)	8 – 23 kg	20/BX
ARJ0120	Video Laryngoscope Sheath (M)	23 – 40 kg	20/BX
ARJ0130	Video Laryngoscope Sheath (L)	>40 kg	20/BX
ARJ13302RU	iS-PF2 Workstation w/wireless connection module		1/EA
ARJ0004RU	Cart		1/EA
ARJ1001RU	Rigid Stylet		1/EA
ARJ020RU	3.5' screen monitor		1/EA
ARJ003RU	Wireless Connection Module		1/EA
ARJ100RU	IS6 Video Laryngoscope Handle		1/EA



teleflex.com/VL

References:

- Crosby ET, Cooper RM, Douglas MJ, et al. *Can J Anaesth.* 1998;45(8):757-776.
- Samsoon GL, Young JR. *Anaesthesia.* 1987;42(5):487-490.
- Lewis SR, Butler AR, Parker J, Cook TM, Smith AF. *Cochrane Database of Systematic Reviews* 2016, Issue 11. Art. No.: CD011136.
- Zhang J, Jiang W, Urdaneta F. *J Comp Eff Res.* 2021 Jul;10(10):831-844. doi: 10.2217/cer-2021-0068. Epub 2021 Apr 27. PMID: 33904779.

Rx Only.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Information in this document is not a substitute for the product Instructions for Use. Refer to package insert for complete warnings, indications, contraindications, precautions, potential complications, and Instructions For Use.

Insighters is a trademark of Shenzhen Insighters Medical Technology Co., Ltd.

Teleflex and the Teleflex logo are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All rights reserved.

© 2023 Teleflex Incorporated. (02/2023) MC-008640