



Rüsch®
TruLite Secure™
 Single-Use Laryngoscope Blade and Handle

Ready When You Are

The Rüsch® TruLite Secure™ Single-Use Laryngoscope from Teleflex is a one-piece solution with a combined blade and handle. This fully disposable, innovative system features reinforced critical components designed to perform in the toughest OR and field conditions. The color-coded, all-in-one blade and handle combination enables quick size selection, while the package design allows light source testing prior to opening.

Bright Light When You Need It Most

The Rüsch® TruLite Secure™ Single-Use Laryngoscope features the latest innovations in blade illumination technology to deliver bright light when it matters most. The LED-in-blade design provides strong, focused, blue-white illumination that is close to the target anatomy for a high definition laryngeal view.

Benefits



The Clinician
 Bright LED light to aid in anatomy visualization in critical situations




Your Institution
 Helps reduce cross-contamination risk vs. reusable devices



The Patient
 A broad range of blade designs and sizes to address diverse patient anatomies

Rüsch®
TruLite Secure™
 Single-Use Laryngoscope Blade and Handle



Color-Coded Sizes
 To ease size identification and selection

LED Light Source
 Provides strong, focused, and reliable illumination

Solid, Metal Construction
 With bevelled blade tip

Universal-Sized Handle
 With textured finger grips for user comfort and ease of use

Battery Compartment
 Contains two alkaline AAA batteries. Safety features help to prevent re-use

Rüsch® TruLite Secure™ Single-Use Laryngoscope Fast Facts

Ideally suited for routine use in the OR or in critical situation settings

Bright blue-white LED light source to support tissue differentiation

Single-use design eliminates workflow and costs associated with reprocessing

One-piece blade and handle system with color-coded sizes for rapid deployment

Rüsch® TruLite Secure™ Single-Use Laryngoscope

PRODUCT CODE	DESCRIPTION	OVERALL LENGTH (MM)	DISTAL WIDTH (MM)	WEIGHT (WITH BATTERIES) (G)	COLOR CODE	CASE QTY
MACINTOSH BLADES						
004671002	Mac 2	100	11.3	108	Turquoise	10
004671003	Mac 3	130	15	125	Yellow	10
004671004	Mac 4	155	15	135	Pink	10
MILLER BLADES						
004670010	Mill 00	67	11.8	97	Beige	10
004670000	Mill 0	78	11.8	98	Purple	10
004670001	Mill 1	103	11.8	104	Orange	10
004670002	Mill 2	155	13.7	125	Grey	10
004670003	Mill 3	195	13.7	137	Green	10
004670004	Mill 4	205	13.7	141	Lavender	10

 teleflex.com/laryngoscopes

Rx Only. Refer to package insert for current warnings, indications, contraindications, precautions, and Instructions for Use.

Teleflex, the Teleflex logo, Rüsch, and TruLite Secure are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. All data current at time of printing (04/2018). Subject to technical changes without further notice. © 2018 Teleflex Incorporated. All rights reserved. MC-000473rev2