




Gel-Bead™ Embolization Spheres

Bioresorbable Gelatin Spheres Designed for Targeted Embolization

Bioresorbable


Gel-Bead™ Embolization Spheres demonstrated degradation beginning at 4 weeks with some spheres still present at 12 weeks.*



4 Weeks
20x magnification, H&E stain

4 Week Histology
Gel-Bead™ 100-300µm Sphere implant site at 4 weeks. As shown above, Gel-Bead™ Spheres began to demonstrate significant degradation (arrows), mediated by the infiltration of inflammatory cells.

*Testing completed within an animal study model by Teleflex. Data on file. Animal study results may or may not be indicative of clinical outcomes in humans.




12 Weeks
20x magnification, H&E stain

12 Week Histology
Gel-Bead™ 100-300µm Sphere implant site at 12 weeks. Tissue surrounding the presumed implant site consisted of fibrous connective tissue containing macrophages and multinucleated giant cells.

Multiple Size Ranges


Available in four size ranges for targeted embolization of hypervascular tumors




Smooth Delivery

Spherical shape designed to provide consistent embolic delivery with predictable sphere size distribution.

Gelatin Spheres



Gelatin Slurry



For illustrative purposes only.

Specifications and Ordering Information

MODEL	SPHERE SIZE	MINIMUM CATHETER I.D.	LABEL COLOR	QTY*
3820	100-300µm	0.020"	Yellow	5
3821	300-500µm	0.020"	Blue	5
3822	500-700µm	0.023"	Red	5
3823	700-1000µm	0.038"	Green	5

*Packaged in quantities of 5 units per box. Individual units are not sold separately.
Gel-Bead™ Embolization Spheres are intended for use in embolization of hypervascular tumors.

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose-driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow®, Deknatel®, Hudson RCI®, LMA®, Pilling®, Rüsch®, UroLift®, and Weck® – trusted brands united by a common sense of purpose.

Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

Regional Offices

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Australia 1300 360 226

Austria +43 (0)1 402 47 72

Belgium +32 (0)2 333 24 60

Canada +1 (0)800 387 9699

China (Shanghai) +86 (0)21 6163 0965

China (Beijing) +86 (0)10 6418 5699

Czech Republic +420 (0)495 759 111

France +33 (0)5 62 18 79 40

Germany +49 (0)711 2090 8000

Greece +30 210 67 77 717

India +91 (0)44 2836 5040

Italy +39 0362 58911

Japan +81 (0)3 6632 3600

Korea +82 2 536 7550

Mexico +52 55 5002 3500

Netherlands +31 (0)88 00 215 00

New Zealand 0800 601 100

Poland +48 22 4624032

Portugal +351 22 541 90 85

Singapore (SEA non-direct sales countries) +65 6439 3000

Slovak Republic +421 (0)3377 254 28

South Africa +27 (0)11 807 4887

Spain +34 918 300 451

Switzerland +41 (0)31 818 40 90

United Kingdom +44 (0)1494 53 27 61

For more information, please visit teleflex.com.

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Teleflex, the Teleflex logo, Arrow, Deknatel, Gel-Bead, Hudson RCI, LMA, Pilling, Rüsch, UroLift, and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries.

Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 02/2020

© 2020 Teleflex Incorporated. All rights reserved. MC-006413