


Axis™ Guidewire

Strength and Support for Stable Access

Heavy-Duty Proximal Shaft
provides distal support and pushability

Smooth Pre-Coating of PTFE
enhances wire delivery

9cm or 12cm Ultra-Floppy Tip
facilitates atraumatic access and device placement



Ordering Information

MODEL	WIRE DIAMETER	WIRE LENGTH	TIP SHAPE	BOX QUANTITY
7370	0.035"	90 cm	Straight, 12cm Floppy	5
7371	0.035"	90 cm	Straight, 9 cm Floppy	5
7372	0.035"	145 cm	Straight, 9 cm Floppy	5
7373	0.035"	260 cm	Straight, 9 cm Floppy	5
7374	0.035"	180 cm	Straight, 9 cm Floppy	5

Packaged in quantities of 10 units per box. Individual units are not sold separately.

The Axis™ guidewire is intended to facilitate the placement of devices during diagnostic and interventional procedures.

VSI Torque Device accessory

MODEL	WIRE DIAMETER	WIRE LENGTH	TIP SHAPE	BOX QUANTITY
8181	0.035"	90 cm	Straight, 12cm Floppy	10

Packaged in quantities of 10 units per box. Individual units are not sold separately.

The Axis™ guidewire is intended to facilitate the placement of devices during diagnostic and interventional procedures.

The VSI Torque Device is designed to be used with guidewires ranging in size from 0.014" to 0.038".

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow®, Deknatel®, Hudson RCI®, LMA®, Pilling®, Rüschi® and Weck® – trusted brands united by a common sense of purpose.

For more information, please visit teleflex.com.

Please see Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Teleflex, the Teleflex logo, Axis, Arrow, Deknatel, Hudson RCI, LMA, Pilling, Rüschi and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All other trademarks and registered trademarks are property of their respective owners.

Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 03/2020.

©2020 Teleflex Incorporated. All rights reserved. MC-006138