

Amplatz SST[™] Guidewire

Extra Stiff for Extra Support

Radiopaque Marker Coil
for enhanced visibility

**Available in Straight,
Floppy and Modified J**
tip configurations

Stainless Steel Core Construction
provides added strength
for increased support

PTFE Coating
for enhanced device tracking

Ordering Information

Amplatz SST™, Straight

MODEL	WIRE DIAMETER	WIRE LENGTH	TIP SHAPE	BOX QUANTITY
7479	0.035"	80 cm	Straight	5
7480	0.035"	150 cm	Straight	5
7481	0.035"	180 cm	Straight	5
7482	0.035"	260 cm	Straight	5

Packaged in quantities of 10 units per box. Individual units are not sold separately.

Amplatz SST™, Modified J

MODEL	WIRE DIAMETER	WIRE LENGTH	TIP SHAPE	BOX QUANTITY
7487	0.035"	150 cm	Modified J	5
7488	0.035"	180 cm	Modified J	5
7489	0.035"	260 cm	Modified J	5

Packaged in quantities of 10 units per box. Individual units are not sold separately.

The Amplatz SST steerable guidewire facilitates placement of a catheter during diagnostic angiography and interventional procedures. The Amplatz SST can be accurately controlled to facilitate navigation through tortuous vessels and/or adjoining side branches.

VSI Torque Device accessory

MODEL	DESCRIPTION	BOX QUANTITY
8181	VSI Torque Device	10

Packaged in quantities of 10 units per box. Individual units are not sold separately.

The VSI Torque Device is intended for use in manipulating a guidewire that is in the vasculature..



Sold Separately

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow®, Deknatel®, Hudson RCI®, LMA®, Pilling®, Rüscher® and Weck® – trusted brands united by a common sense of purpose.

For more information, please visit teleflex.com.

Please see Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Teleflex, the Teleflex logo, and AmplatzSST are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All other trademarks and registered trademarks are property of their respective owners.

Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 04/2020.

©2020 Teleflex Incorporated. All rights reserved. MC-006152