

VSI StraitSet™  
Micro-Introducer

GrebSet®  
Micro-Introducer

VSI StraitSet™ Micro-Introducer Kit  
GrebSet® Micro-Introducer Kit


# GrebSet<sup>®</sup>

micro-introducer kit

## A Unique Angle on Micro Access


### 5 Fr. Coaxial Design

- Angled distal tip for steerability and control of 0.038" guidewire
- Hydrophilic coating for smooth delivery
- Embedded radiopaque marker band for visibility
- Locking trocar needle option to help prevent back-out during insertion




#### Full Length Stainless Steel Braid

provides excellent flexibility, pushability and kink resistance



Radiopaque Marker



## GrebSet<sup>®</sup> Micro-Introducer Kit Components

- 5 Fr. braided sheath with hydrophilic coating available in 30 cm and 45 cm lengths
- Stiffen dilator
- .018" nitinol mandrel guidewire with angled, tungsten coil tip available in 70 cm and 100 cm lengths

### Needle Options

- 21 Ga. x 15 cm marked trocar needle with echogenic tip and depth stop
- 21 Ga. x 7 cm echogenic needle
- 21 Ga. x 15 cm marked trocar needle with echogenic tip and locking stylet

# VSI StraitSet<sup>™</sup>

micro-introducer kit

## Straight Simplicity

### 6 Fr. Triaxial Design

- Straight sheath for simple guidewire introduction, allows dual placement of a 0.038" and 0.018" guidewire
- Hydrophilic or siliconized coating for smooth delivery
- Embedded radiopaque marker band for visibility
- Locking trocar needle option to help prevent back-out during insertion



Radiopaque Marker



## VSI StraitSet<sup>™</sup> Micro-Introducer Kit Components

- 6 Fr. x 20 cm sheath with hydrophilic coating or silicone wipe
- Dilator
- Metal stiffening cannula
- .018" x 60 cm nitinol or stainless steel mandrel guidewire with tungsten tip

### 0.038" wire option:

- 0.038" x 145 cm PTFE 3 mm "J"/straight guidewire

### Needle options

- 21 Ga. x 15 cm marked trocar needle with echogenic tip and locking stylet
- 21 Ga. x 7 cm echogenic needle

# Ordering Information

## GrebSet® Micro-Introducer Kit

MODEL	INTRODUCER SIZE	SHEATH LENGTH	NEEDLE	SHEATH COATING	WIRE DIAMETER	WIRE LENGTH	WIRE MANDREL	WIRE COIL TIP	BOX QUANTITY
7941	5 Fr.	45 cm	21G x 15 cm Echogenic Trocar Needle with Depth Markings and	Hydrophilic	0.018"	100 cm	Nitinol	Tungsten	10
7950	5 Fr.	30 cm	21G x 7 cm Echogenic Needle	Hydrophilic	0.018"	70 cm	Nitinol	Tungsten	10
7951	5 Fr.	30 cm	21G x 15 cm Marked, Echogenic Trocar Needle	Hydrophilic	0.018"	70 cm	Nitinol	Tungsten	10
7952	5 Fr.	30 cm	21G x 15 cm Echogenic Trocar Needle with Depth Markings and Locking Stylet	Hydrophilic	0.018"	70 cm	Nitinol	Tungsten	10

Packaged in quantities of 10 units per box. Individual units are not sold separately.

The GrebSet® Micro-Introducer Kit is indicated for use in percutaneous introduction and placement of catheters and guidewires.

## VSI StraitSet™ Micro-Introducer Kit

MODEL	SHEATH O.D.	SHEATH I.D.	SHEATH LENGTH	NEEDLE	SHEATH COATING	0.018" WIRE LENGTH	0.018" WIRE MANDREL	0.018" WIRE COIL TIP	0.018" WIRE	BOX QUANTITY
7960	6 Fr.	4 Fr.	20 cm	21G x 7 cm Echogenic Needle	Hydrophilic	60 cm	Nitinol	Tungsten	N/A	10
7961	6 Fr.	4 Fr.	20 cm	21G x 15 cm Echogenic Trocar Needle with Depth Markings and Locking Stylet	Hydrophilic	60 cm	Nitinol	Tungsten	N/A	10
7962	6 Fr.	4 Fr.	20 cm	21G x 15 cm Echogenic Trocar Needle with Depth Markings and Locking Stylet	Hydrophilic	60 cm	Stainless Steel	Tungsten	N/A	10
7964	6 Fr.	4 Fr.	20 cm	21G x 15 cm Echogenic Trocar Needle with Depth Markings and Locking Stylet	Hydrophilic	60 cm	Nitinol	Tungsten	YES	10
7965	6 Fr.	4 Fr.	20 cm	21G x 15 cm Echogenic Trocar Needle with Depth Markings and Locking Stylet	Hydrophilic	60 cm	Stainless Steel	Tungsten	YES	10
7971	6 Fr.	4 Fr.	20 cm	21G x 15 cm Echogenic Trocar Needle with Depth Markings and Locking Stylet	Siliconized	60 cm	Nitinol	Tungsten	N/A	10

Packaged in quantities of 10 units per box. Individual units are not sold separately.

The VSI StraitSet™ Micro-Introducer Kit is indicated for use in percutaneous introduction and placement of catheters and guidewires.

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Teleflex, the Teleflex logo, GrebSet, and VSI StraitSet are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries.

Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 05/2020

© 2020 Teleflex Incorporated. All rights reserved. MC-006343