



ACCESSORY DEVICES

for Your Vein Practice

CHOICE • VALUE • CONVENIENCE



YOUR PREFERENCES FOR YOUR PROCEDURES

VSI is known for working closely with physicians to design products based on their specific clinical needs, and offers a wide breadth of accessory devices for venous ablation procedures.

7F INTRODUCER KITS

- 7F sheath available in 7cm, 11cm or 12cm lengths
- 21G or 20G, 7cm needle and 0.018" guidewire with nitinol mandrel and stainless steel tip



NON-COATED SHEATHS



COATED SHEATHS

NON-COATED SHEATHS

Model	Introducer Size	Sheath Length	Needle Tip	Needle Length	Sheath Coating	Guidewire Diameter	Guidewire Length	Guidewire Mandrel	Guidewire Coil Tip	Box Quantity
8981	7F	7cm	20G Smooth	7cm	Non-coated	N/A	N/A	N/A	N/A	5
8982	7F	11cm	20G Smooth	7cm	Non-coated	N/A	N/A	N/A	N/A	5
8983*	7F	7cm	20G Smooth	7cm	Non-coated	0.018"	40cm	Nitinol	Stainless Steel	5
8984*	7F	11cm	20G Smooth	7cm	Non-coated	0.018"	40cm	Nitinol	Stainless Steel	5

Individual units are not sold separately. *Guidewire individually pouched.

7F Introducer Kits are intended for use to facilitate the introduction of guidewires, catheters and other accessory medical devices through the skin into a vein or artery and minimize blood loss associated with such introduction.

COATED SHEATHS

Model	Introducer Size	Sheath Length	Needle Tip	Needle Length	Sheath Coating	Guidewire Diameter	Guidewire Length	Guidewire Mandrel	Guidewire Coil Tip	Box Quantity
8977	7F	12cm	21G Smooth	7cm	Siliconized	0.018"	40cm	Nitinol	Stainless Steel	5
8979	7F	7cm	21G Smooth	7cm	Siliconized	0.018"	40cm	Nitinol	Stainless Steel	5

Individual units are not sold separately.

The VSI Micro-HV Introducer Kits are intended for use in percutaneous procedures to introduce catheters and other intravascular devices into the vasculature.

VSI STANDARD PROCEDURE PACKS

Standard procedure packs contain pre-determined, frequently requested components for both endovenous laser ablation and radiofrequency ablation procedures.



ECONOMY PROCEDURE PACK WITH CHLORAPREP® - Model 7785

Syringes, Needles and Scalpels

- 1 30ml luer lock syringe
- 1 5ml luer lock syringe
- 1 25 gauge x 1½" needle
- 1 22 gauge x 3½" spinal needle
- 1 Disposable #11 blade scalpel with safety guard

Accessory Items

- 1 Lidocaine HCL USP 1% – 10mg/ml, 20ml
- 4 Ultrasound gel – 20g each
- 1 ChloraPrep applicator 10.5ml
- 1 250cc solution bowl
- 2 Latex free rubberbands
- 1 Tray 14" x 8" x 10"

Gowns, Drapes, Towels and Dressings

- 1 Prep table cover
- 1 XL disposable gown
- 4 Utility drapes with tape - 15" x 26"
- 1 ½ drape – 60" x 44"
- 1 Drape – 60" x 76"
- 10 4" x 4" gauze pads
- 4 OR towels
- 1 Ultrasound probe cover - 6" x 96"
- 1 Bandage – 4" x 5yds

BASIC PROCEDURE PACK - Model 7790

Syringes, Needles and Scalpels

- 1 30ml luer lock syringe
- 1 10ml luer lock syringe
- 1 5ml luer lock syringe
- 1 25 gauge x 1½" needle
- 1 19 gauge x 1½" needle
- 1 Disposable #11 blade scalpel with safety guard

Accessory Items

- 1 Lidocaine HCL USP 1% – 10mg/ml, 20ml
- 3 Ultrasound gel – 20g each
- 1 Iodine prep solution – 2oz
- 1 Prep stick sponge
- 2 250cc solution bowls
- 1 Skin marking pen with 6" flexible ruler
- 2 Rubberbands #32
- 1 Tray 11 x 15 procedure set-up

Gowns, Drapes, Towels and Dressings

- 1 Prep table cover
- 2 XL disposable gowns
- 1 ½ drape – 40" x 58"
- 1 Split drape with tape – 77" x 120"
- 10 4" x 4" gauze pads
- 4 OR towels
- 1 Ultrasound probe cover

COMPREHENSIVE PROCEDURE PACK - Model 7791

Syringes, Needles and Scalpels

- 2 30ml luer lock syringes
- 2 10ml luer lock syringes
- 1 5ml luer lock syringe
- 1 25 gauge x 1½" needle
- 1 20 gauge x 3½" needle
- 1 19 gauge x 1½" needle
- 1 Disposable #11 blade scalpel with safety guard
- 1 Sharps holder (Short Stop)

Accessory Items

- 1 Lidocaine HCL USP 1% – 10mg/ml, 20ml
- 4 Ultrasound gel – 20g each
- 1 Iodine prep solution – 2oz
- 1 Prep stick sponge
- 1 500cc utility bowl
- 1 250cc solution bowl
- 1 Skin marking pen with 6" flexible ruler
- 2 Rubberbands #32
- 1 Tray 11 x 15 procedure set-up

Gowns, Drapes, Towels and Dressings

- 1 Prep table cover
- 2 XL disposable gowns
- 1 ½ drape – 40" x 58"
- 1 Split drape with tape – 77" x 120"
- 4 Utility drapes – 15" x 26"
- 20 4" x 4" gauze pads
- 4 OR towels
- 1 Ultrasound probe cover
- 1 Steri strip – ½" x 4"
- 1 Self-adherent wrap – 4" x 5yds

COMPREHENSIVE PROCEDURE PACK WITH CHLORAPREP® - Model 7792

Syringes, Needles and Scalpels

- 2 30ml luer lock syringes
- 2 10ml luer lock syringes
- 1 5ml luer lock syringe
- 1 25 gauge x 1½" needle
- 1 20 gauge x 3½" needle
- 1 19 gauge x 1½" needle
- 1 Disposable #11 blade scalpel with safety guard
- 1 Sharps holder (Short Stop)

Accessory Items

- 1 Lidocaine HCL USP 1% – 10mg/ml, 20ml
- 4 Ultrasound gel – 20g each
- 1 ChloraPrep applicator 10.5ml
- 1 250cc solution bowl
- 1 Skin marking pen with 6" flexible ruler
- 2 Rubberbands #32
- 1 Tray 11 x 15 procedure set-up

Gowns, Drapes, Towels and Dressings

- 1 Prep table cover
- 2 XL disposable gowns
- 1 ½ drape – 40" x 58"
- 1 Split drape with tape – 77" x 120"
- 4 Utility drapes – 15" x 26"
- 20 4" x 4" gauze pads
- 4 OR towels
- 1 Ultrasound probe cover
- 1 Steri strip ½" x 4"
- 1 Self-adherent wrap – 4" x 5yds

Model	Description
7785	Economy Procedure Pack with ChloraPrep®
7790	Basic Procedure Pack
7791	Comprehensive Procedure Pack
7792	Comprehensive Procedure Pack with ChloraPrep®

VSI RF ABLATION PROCEDURE PACKS AND KITS



RADIOFREQUENCY ABLATION PROCEDURE PACK - Model 7795

Radiofrequency Ablation Procedure Pack 7F 7cm Micro-HV Kit - Model 7796

Radiofrequency Ablation Procedure Pack 7F 12cm Micro-HV Kit - Model 7797

Radiofrequency Ablation Procedure Pack 7F 7cm Introducer Kit - Model 7788

Radiofrequency Ablation Procedure Pack 7F 11cm Introducer Kit - Model 7789

Syringes, Needles and Scalpels

- 2 30ml luer lock syringes
- 1 20ml luer lock syringe
- 1 10ml luer lock syringe
- 1 27 gauge x 1 ¼" needle
- 1 20 gauge x 3 ½" spinal needle
- 1 18 gauge x 1 ½" needle
- 1 Disposable #11 blade scalpel with safety guard

Accessory Items

- 3 Ultrasound gel – 20ml
- 2 Prep stick sponges
- 1 Tray – 9" x 5" x 3"
- 1 500cc utility bowl
- 1 250cc solution bowl
- 1 Skin marking pen
- 2 Latex free rubber bands
- 1 Zip lock bag

Gowns, Drapes, Towels and Dressings

- 1 Prep table cover
- 2 XL disposable gowns
- 2 Utility drapes with tape - 15" x 26"
- 1 ½ drape – 60" x 44"
- 1 U-drape – 80" x 124" – 4" x 44"
- 30 4" x 4" gauze pads
- 3 OR towels, blue
- 2 ABS towels, white
- 1 Towel clamp
- 1 Ultrasound probe cover – 6" x 96"
- 1 Steri strip ½" x 4"

BASIC RADIOFREQUENCY ABLATION PROCEDURE PACK - Model 7798

Basic Radiofrequency Ablation Procedure Pack 7F 7cm Introducer Kit - Model 7786

Basic Radiofrequency Ablation Procedure Pack 7F 11cm Introducer Kit - Model 7787

Syringes, Needles and Scalpels

- 1 30ml luer lock syringe
- 1 20ml luer lock syringe
- 1 10ml luer lock syringe
- 1 27 gauge x 1 ¼" needle
- 1 20 gauge x 3 ½" spinal needle
- 1 18 gauge x 1 ½" needle
- 1 Disposable #11 blade scalpel with safety guard

Accessory Items

- 2 Ultrasound gel – 20ml
- 2 Prep stick sponges
- 1 Tray – 9"x5"x2"
- 1 500cc utility bowl
- 1 250cc utility bowl
- 2 Latex free rubber bands
- 1 Zip lock bag

Gowns, Drapes, Towels and Dressings

- 1 Prep table cover
- 2 XL disposable gowns
- 1 U-Drape - 80" x 124" – 4" x 44"
- 30 4"x4" gauze pads
- 2 OR towels, blue
- 1 ABS towels, white
- 1 Towel clamp
- 1 Ultrasound probe cover – 6" x 96"

Model	Description
7795	Radiofrequency Ablation Procedure Pack
7798	Basic Radiofrequency Ablation Procedure Pack
Radiofrequency Ablation Procedure Pack Kits	
7796	Radiofrequency Ablation Procedure Pack with 7F 7cm Micro-HV Kit (Coated)
7797	Radiofrequency Ablation Procedure Pack with 7F 12cm Micro-HV Kit (Coated)
7788	Radiofrequency Ablation Procedure Pack 7F 7cm Introducer Kit (Non-Coated)
7789	Radiofrequency Ablation Procedure Pack 7F 11cm Introducer Kit (Non-Coated)
Basic Radiofrequency Ablation Procedure Pack Kits	
7786	Basic Radiofrequency Ablation Procedure Pack 7F 7cm Introducer Kit (Non-Coated)
7787	Basic Radiofrequency Ablation Procedure Pack 7F 11cm Introducer Kit (Non-Coated)

VSI SELECT NEEDLE™

This 18G x 7cm echogenic needle promotes highly visible access under ultrasound.



Model	Description	Gauge	Length	Tip	Box Quantity
7650	VSI Select Needle	18G	7cm	Echogenic	10

Individual units are not sold separately.

The VSI Select Needle is intended for use as a percutaneous entry needle.

0.018" AND 0.025" GUIDEWIRES

0.018" and 0.025" guidewire configurations are available in a variety of length, mandrel, and tip choices. 0.018" guidewires promote venous access while 0.025" guidewires facilitate placement.



Model	Diameter	Length	Mandrel	Tip
7311	0.018"	100cm	Nitinol	Stainless Steel
7312	0.018"	120cm	Nitinol	Stainless Steel
7313	0.018"	150cm	Nitinol	Stainless Steel

Each component is individually pouched and sold individually.

Model	Diameter	Length	Tip	Box Quantity
7385	0.025"	150cm	Double-ended, Straight and 1.5mm J-tip	5

Individual units are not sold separately.

The VSI 0.018" Guidewire is indicated for percutaneous entry of peripheral vessels using the Seldinger Technique. This device is not intended for use in the coronary or cerebral vasculature. The VSI .025" guidewire is intended to facilitate the placement of devices during diagnostic and interventional procedures.

AUTO-FILL® SYRINGE KIT

The Auto-Fill syringe kit offers a simple fluid delivery system for local anesthesia. A polycarbonate control syringe, IV extension line and three-way check valve ensure maximum efficiency and allow the delivery of multiple fluid doses without the need to manually reload the syringe.



Specifications:
• 12cc polycarbonate syringe (accuracy ±5%)
• IV tubing (108", 274cm) with spike, three-way check valve, and roller clamp
• 36" (91.44cm) extension tubing

Model	Description	Quantity
7600	Auto-Fill Syringe Kit	10

Individual units are not sold separately.

The Auto-Fill Syringe Kit is indicated for introduction of dilute lidocaine solutions into the subcutaneous tissue for the purposes of local anesthesia.

KLEIN PUMP® TUBING

The Klein Pump Tubing is a disposable small bore tubing for infusion of anesthetic with an infusion pump.



Model	Description	Quantity
CVT-10VS	Klein CV Infiltration Set - Single Spike*	10

Individual units are not sold separately.

The Klein Pump is indicated to cause a flow of fluid from an IV bag into a patient in a manner controlled manually by a healthcare professional. The Klein Pump is not intended to be used as an IV infusion pump.

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Auto-Fill is a registered trademark of Vascular Solutions, Inc. ChloroPrep is a registered trademark of CareFusion 2200, Inc. VSI Select Needle is a trademark of Vascular Solutions, Inc. Klein Pump is a registered trademark of HK Surgical. Ziploc is a registered trademark of S.C. Johnson & Son, Inc.

©2017 Vascular Solutions, Inc. All rights reserved. ML2963 Rev. C 06/17



Vascular Solutions, Inc.

6464 Sycamore Court North
Minneapolis, Minnesota 55369 USA

Customer Service: 888.240.6001

customerservice@vasc.com

www.vasc.com