




50SC06L01 Pilling™ Microlaryngoscopy Solid Container System


Pilling™

## ENT Instruments

Instrument Configurations for Microlaryngoscopy,  
Bronchoscopy, Esophagoscopy, and Panendoscopy

# Microlaryngoscopy Configurations

## Solid Container




- Solid bottom is compatible with Steam, EtO, and Sterizone VP4 sterilization methods
- Utilizes round sterilization filters

### Solid Containers

Three basket levels hold 6 laryngoscopes, plus instruments and light cables.

|                          |           |                                                                                                                                                                                                                                                                                   |
|--------------------------|-----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> | 50SC06L01 | Pilling™ Micro Laryngoscopy Container, 27" x 11" x 6", with Starter Kit<br>Three Basket Levels: bottom 3" scope basket, middle & upper 1" instrument baskets<br>Starter Kit includes: 40 disposable filters, 20 tamper evident seals, 20 load cards, and 20 load indicator strips |
|--------------------------|-----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## Wrapped Container



- Features two tiers: bottom tier holds six laryngoscopes, top tier stores instrumentation
- Perforated container requires wrapping for sterilization

### Wrapped Container Options

Two-tier system. Bottom level holds 6 laryngoscopes and light cables. Top tray holds instrumentation.

|                          |           |                                                                         |
|--------------------------|-----------|-------------------------------------------------------------------------|
| <input type="checkbox"/> | P88097    | Double-Decker ENT Tray w/Inserts and Standard Mats, 23" x 11" x 5"      |
| <input type="checkbox"/> | P88044    | Double-Decker ENT Tray w/Inserts Only, 23" x 11" x 5"                   |
| <input type="checkbox"/> | 087396 x2 | Ultrasoft Silicone Spiked Mat, 22" L x 11" W – to be paired with P88044 |



## Laryngoscopes

Light carrier included.

|                          |        |                                                                      |
|--------------------------|--------|----------------------------------------------------------------------|
| <input type="checkbox"/> | 532215 | <b>NEW</b> Lindholm-style Laser Laryngoscope, 38 mm x 150 mm         |
| <input type="checkbox"/> | 523217 | <b>NEW</b> Lindholm-style Laser Laryngoscope, 38 mm x 170 mm         |
| <input type="checkbox"/> | 532221 | <b>NEW</b> Adult Dedo Laser Micro Laryngoscope, 19 mm x 178 mm       |
| <input type="checkbox"/> | 532206 | <b>NEW</b> Jako Type 3 Laryngoscope, 29 mm x 178 mm                  |
| <input type="checkbox"/> | 532196 | <b>NEW</b> Garrett-Ossoff-Pilling Laryngoscope, 25 mm x 188 mm       |
| <input type="checkbox"/> | 532191 | <b>NEW</b> Ossoff-Pilling Laser Laryngoscope, 25 mm x 168mm          |
| <input type="checkbox"/> | 532190 | <b>NEW</b> Ossoff-Pilling Laser Laryngoscope, 19 mm x 168 mm         |
| <input type="checkbox"/> | 532010 | <b>NEW</b> Holinger Anterior Commissure Laryngoscope, 10 mm x 170 mm |

## Laryngeal Instruments

### Cup Forceps – Standard (22 cm – 28 cm)

|                          |        |                                      |
|--------------------------|--------|--------------------------------------|
| <input type="checkbox"/> | 506455 | Jako 2 mm Biopsy Straight, 22 cm     |
| <input type="checkbox"/> | 506460 | Jako 2 mm Biopsy Angled Left, 22 cm  |
| <input type="checkbox"/> | 506461 | Jako 2 mm Biopsy Angled Right, 22 cm |
| <input type="checkbox"/> | 506462 | Jako 2 mm Angled Up, 22 cm           |
| <input type="checkbox"/> | 505100 | Jackson 4 mm Straight, 28 cm         |
| <input type="checkbox"/> | 505104 | Jackson 4 mm Angled Up, 28 cm        |
| <input type="checkbox"/> | 505102 | Jackson 6 mm Straight, 28 cm         |
| <input type="checkbox"/> | 505106 | Jackson 6 mm Angled Up, 28 cm        |

### Micro-Cup Forceps (23 cm)

|                          |        |                      |
|--------------------------|--------|----------------------|
| <input type="checkbox"/> | 506862 | Phono 1 mm Straight  |
| <input type="checkbox"/> | 506864 | Phono 1 mm Angled Up |
| <input type="checkbox"/> | P28526 | Phono 1 mm Right     |
| <input type="checkbox"/> | P28527 | Phono 1 mm Left      |

### Grasping Forceps – Standard (22 cm)

|                          |        |                                |
|--------------------------|--------|--------------------------------|
| <input type="checkbox"/> | 506470 | Jako Alligator Curved Left     |
| <input type="checkbox"/> | 506471 | Jako Alligator Curved Right    |
| <input type="checkbox"/> | 506472 | Jako Alligator Straight        |
| <input type="checkbox"/> | 506950 | Laryngeal Triangular Jaw Right |
| <input type="checkbox"/> | 506952 | Laryngeal Triangular Jaw Left  |

### Micro-Grasping Forceps (23 cm)

|                          |        |                                       |
|--------------------------|--------|---------------------------------------|
| <input type="checkbox"/> | 506856 | Phono 1 mm Double Action with Suction |
| <input type="checkbox"/> | 506858 | Phono 1 mm Angled Up                  |
| <input type="checkbox"/> | 506860 | Phono 1 mm Straight                   |
| <input type="checkbox"/> | P28531 | Phono 1 mm Right                      |
| <input type="checkbox"/> | P28532 | Phono 1 mm Left                       |

### Scissors – Standard (22 cm)

|                          |        |                          |
|--------------------------|--------|--------------------------|
| <input type="checkbox"/> | 506475 | Jako Curved Left         |
| <input type="checkbox"/> | 506476 | Jako Curved Right        |
| <input type="checkbox"/> | 506477 | Jako Slight Up Angle     |
| <input type="checkbox"/> | 506478 | Jako Straight Vertical   |
| <input type="checkbox"/> | 506479 | Jako Straight Horizontal |

### Suction Tubes (20 cm – 35 cm)

|                          |        |                                                                       |
|--------------------------|--------|-----------------------------------------------------------------------|
| <input type="checkbox"/> | 507010 | 12 Fr Jackson Open End, 20 cm                                         |
| <input type="checkbox"/> | 507030 | 12 Fr Jackson Velvet Eye, 20 cm                                       |
| <input type="checkbox"/> | 507065 | 12 Fr Holinger Velvet Eye, 35 cm                                      |
| <input type="checkbox"/> | 507060 | 9 Fr Holinger Velvet Eye, 20 cm                                       |
| <input type="checkbox"/> | 507122 | Jako 7 Fr with Thumb Control                                          |
| <input type="checkbox"/> | 507121 | Jako Suction Tube 6 Fr with Thumb Control                             |
| <input type="checkbox"/> | 507120 | Jako 5 Fr with Thumb Control                                          |
| <input type="checkbox"/> | 507045 | Insulated 12 Fr Suction Cautery Tube, 22 cm w/10 ft Bovie Cord & Plug |

### Knives, Probes, and Elevators (27 cm)

|                          |        |                                                                     |
|--------------------------|--------|---------------------------------------------------------------------|
| <input type="checkbox"/> | 506834 | <b>NEW</b> Ossoff Micro Flap Elevator, 30°, 1.0 mm                  |
| <input type="checkbox"/> | 506835 | <b>NEW</b> Ossoff Micro Flap Elevator, 30°, 1.5 mm                  |
| <input type="checkbox"/> | 506836 | <b>NEW</b> Ossoff Micro Flap Elevator, 30°, 2.0 mm                  |
| <input type="checkbox"/> | 506837 | <b>NEW</b> Ossoff Micro Flap Elevator Set, One of Each Size         |
| <input type="checkbox"/> | 506830 | Phono Elevator 30° Blunt, 3.0 mm                                    |
| <input type="checkbox"/> | 506832 | Phono Elevator 60° Blunt, 3.0 mm                                    |
| <input type="checkbox"/> | 506673 | Jako Blunt Probe Angled                                             |
| <input type="checkbox"/> | 506669 | Jako Blunt Probe Straight                                           |
| <input type="checkbox"/> | 506880 | Phono Sickle Knife Straight                                         |
| <input type="checkbox"/> | 506882 | Phono Sickle Knife Curved                                           |
| <input type="checkbox"/> | 506888 | Universal Ring Handle with Camlock – for use with knives and probes |

### Ventilation

|                          |        |                                                              |
|--------------------------|--------|--------------------------------------------------------------|
| <input type="checkbox"/> | 502600 | Venturi Ventilating Needle, Blunt 14 Gauge, Right Angle, 4"  |
| <input type="checkbox"/> | 502594 | Venturi Ventilating Needle, Blunt 14 Gauge, Straight, 5-3/8" |
| <input type="checkbox"/> | 522240 | Clip for Venturi Ventilating Needles                         |

## Accessories

### Light Cable

|                          |         |                                                   |
|--------------------------|---------|---------------------------------------------------|
| <input type="checkbox"/> | P358PRA | Pilling™ to Pilling™ Blue Fiber Optic Cable, 8 ft |
| <input type="checkbox"/> | A358PRA | ACMI to Pilling™ Blue Fiber Optic Cable, 8 ft     |

### Solid Container Accessories

|                          |          |                                   |
|--------------------------|----------|-----------------------------------|
| <input type="checkbox"/> | 50SCS01B | Blue Seal for Steam-EtO, 1000pk   |
| <input type="checkbox"/> | 50SCL01  | Load Card for Steam-EtO, 1000pk   |
| <input type="checkbox"/> | 50SCF01  | Filters, Cellulosic, 7.5", 1000pk |

### Teeth Guards

|                          |        |                                     |
|--------------------------|--------|-------------------------------------|
| <input type="checkbox"/> | 502275 | Disposable Teeth Guards 24 Pair/Box |
|--------------------------|--------|-------------------------------------|



To view Laryngoscopy products go to:  
**myTeleflex.com**

# Bronchoscopy Checklist



Bronchoscopy Configuration in Tray P37331

To view Bronchoscopy products go to:  
[myTeleflex.com](http://myTeleflex.com)



## Bronchoscopes

### Holinger-Type Laser Bronchoscopes

*Light carrier and window plug included.*

*Contain ports for flexible instruments and suction.*

|                          |         |                                           |
|--------------------------|---------|-------------------------------------------|
| <input type="checkbox"/> | 522556  | Improved Holinger, 4 mm x 25 cm           |
| <input type="checkbox"/> | 523120P | Pilling™ Laser Bronchoscope, 5 mm x 30 cm |
| <input type="checkbox"/> | 523130P | Pilling™ Laser Bronchoscope, 6 mm x 35 cm |
| <input type="checkbox"/> | 523140P | Pilling™ Laser Bronchoscope, 7 mm x 40 cm |
| <input type="checkbox"/> | 523150P | Pilling™ Laser Bronchoscope, 8 mm x 40 cm |
| <input type="checkbox"/> | 523160P | Pilling™ Laser Bronchoscope, 9 mm x 40 cm |

### Alternative: Holinger Bronchoscopes

*Light carrier and window plug included.*

|                          |        |                                     |
|--------------------------|--------|-------------------------------------|
| <input type="checkbox"/> | 522578 | Holinger Bronchoscope, 4 mm x 30 cm |
| <input type="checkbox"/> | 522579 | Holinger Bronchoscope, 5 mm x 30 cm |
| <input type="checkbox"/> | 522580 | Holinger Bronchoscope, 5 mm x 35 cm |
| <input type="checkbox"/> | 522581 | Holinger Bronchoscope, 5 mm x 40 cm |
| <input type="checkbox"/> | 522582 | Holinger Bronchoscope, 5 mm x 45 cm |
| <input type="checkbox"/> | 522583 | Holinger Bronchoscope, 6 mm x 35 cm |
| <input type="checkbox"/> | 522584 | Holinger Bronchoscope, 6 mm x 40 cm |
| <input type="checkbox"/> | 522585 | Holinger Bronchoscope, 7 mm x 40 cm |
| <input type="checkbox"/> | 522586 | Holinger Bronchoscope, 8 mm x 40 cm |
| <input type="checkbox"/> | 522587 | Holinger Bronchoscope, 9 mm x 40 cm |

## Accessories

### Light Cable

|                          |         |                                                   |
|--------------------------|---------|---------------------------------------------------|
| <input type="checkbox"/> | P358PRA | Pilling™ to Pilling™ Blue Fiber Optic Cable, 8 ft |
| <input type="checkbox"/> | A358PRA | ACMI to Pilling™ Blue Fiber Optic Cable, 8 ft     |

### Visualization


|                          |        |                              |
|--------------------------|--------|------------------------------|
| <input type="checkbox"/> | 502540 | Fluvog Slide Assembly        |
| <input type="checkbox"/> | 524050 | Telescope 0°, 4.0 mm X 45 cm |
| <input type="checkbox"/> | 524022 | Telescope 0°, 2.7 mm X 34 cm |
| <input type="checkbox"/> | 524037 | 4.0 mm Telescope Guide       |
| <input type="checkbox"/> | 524036 | 2.7 mm Telescope Guide       |

## Tray Setup

|                          |           |                                                     |
|--------------------------|-----------|-----------------------------------------------------|
| <input type="checkbox"/> | P37331    | Single-Tier ENT Tray w/o Brackets, 11" x 23" x 3.5" |
| <input type="checkbox"/> | P37306 x3 | Universal Keyhole Bracket 15 Position, 10"          |
| <input type="checkbox"/> | P37305 x2 | Accessory Box w/Hardware, 1.5" x 3.5" x 1"          |
| <input type="checkbox"/> | 087396    | Ultrasoft Silicone Spiked Mat, 11" W x 22" L        |



# Esophagoscopy Checklist



Esophagoscopy Configuration in Tray P37331

## Esophagoscopes

### Adult Esophagoscopes – Oval

*Light carrier included.*

|                          |        |                                                         |
|--------------------------|--------|---------------------------------------------------------|
| <input type="checkbox"/> | 523640 | Pillsbury Jesberg Esophagoscope, 6 mm x 8 mm x 32 cm    |
| <input type="checkbox"/> | 523642 | Pillsbury Jesberg Esophagoscope, 6.5 mm x 10 mm x 35 cm |
| <input type="checkbox"/> | 523644 | Pillsbury Jesberg Esophagoscope, 7 mm x 10 mm x 42 cm   |
| <input type="checkbox"/> | 523646 | Pillsbury Jesberg Esophagoscope, 8 mm x 12 mm x 45 cm   |
| <input type="checkbox"/> | 523648 | Pillsbury Jesberg Esophagoscope, 10 mm x 14 mm x 29 cm  |
| <input type="checkbox"/> | 523649 | Pillsbury Jesberg Esophagoscope, 10 mm x 14 mm x 50 cm  |
| <input type="checkbox"/> | 523677 | Forbes Esophageal Specula Adult, 22.2 cm                |

### Alternative – Holinger Esophagoscopes – Round

*Light carrier included.*

|                          |        |                                      |
|--------------------------|--------|--------------------------------------|
| <input type="checkbox"/> | 523571 | Holinger Esophagoscope, 4 mm x 30 cm |
| <input type="checkbox"/> | 523572 | Holinger Esophagoscope, 6 mm x 30 cm |
| <input type="checkbox"/> | 523573 | Holinger Esophagoscope, 7 mm x 30 cm |
| <input type="checkbox"/> | 523574 | Holinger Esophagoscope, 9 mm x 30 cm |
| <input type="checkbox"/> | 523575 | Holinger Esophagoscope, 5 mm x 35 cm |
| <input type="checkbox"/> | 523577 | Holinger Esophagoscope, 7 mm x 45 cm |
| <input type="checkbox"/> | 523579 | Holinger Esophagoscope, 9 mm x 53 cm |
| <input type="checkbox"/> | 523581 | Holinger Esophagoscope, 8 mm x 45 cm |
| <input type="checkbox"/> | 523582 | Holinger Esophagoscope, 9 mm x 45 cm |

## Accessories

### Light Cable

|                          |         |                                                   |
|--------------------------|---------|---------------------------------------------------|
| <input type="checkbox"/> | P358PRA | Pilling™ to Pilling™ Blue Fiber Optic Cable, 8 ft |
| <input type="checkbox"/> | A358PRA | ACMI to Pilling™ Blue Fiber Optic Cable, 8 ft     |


### Tray Setup

|                          |           |                                                     |
|--------------------------|-----------|-----------------------------------------------------|
| <input type="checkbox"/> | P37331    | Single-Tier ENT Tray w/o Brackets, 11" x 23" x 3.5" |
| <input type="checkbox"/> | P37306 x2 | Universal Keyhole Bracket 15 Position, 10"          |
| <input type="checkbox"/> | 087396    | Ultrasoft Silicone Spiked Mat, 11" W x 22" L        |




To view Esophagoscopy products go to:  
**myTeleflex.com**

# Bronchoscopy and Esophagoscopy Instrument Checklists



Bronchoscopy & Esophagoscopy  
Instrument Configuration in  
Tray P37304



Bronchoscopy & Esophagoscopy  
XL Instrument Configuration in  
Tray P37319

## Standard-Length Instruments (up to 50 cm)

### Flexible Instruments for Bronchoscopes

Primarily for side-channel use with Improved Holinger Bronchoscopes.

|                          |        |                                                   |
|--------------------------|--------|---------------------------------------------------|
| <input type="checkbox"/> | 505670 | Flexible Grasper Forceps, 1 mm                    |
| <input type="checkbox"/> | 505674 | Flexible Biopsy Cup Forceps, 1 mm                 |
| <input type="checkbox"/> | 505680 | Flexilble Alligator Grasper Forceps, 2 mm x 40 cm |
| <input type="checkbox"/> | 505686 | Flexible Alligator Grasper Forceps, 2 mm x 50 cm  |
| <input type="checkbox"/> | 505690 | Flexible Circular Cup Forceps, 2 mm x 50 cm       |

### Cup Biopsy Forceps

|                          |        |                                                  |
|--------------------------|--------|--------------------------------------------------|
| <input type="checkbox"/> | 505654 | Pediatric Straight Cup Forceps, 35 cm            |
| <input type="checkbox"/> | 505132 | Roberts Straight 4 mm Cup Forceps, 50 cm         |
| <input type="checkbox"/> | 505144 | Roberts Angled 4 mm Cup Forceps, 50 cm           |
| <input type="checkbox"/> | 505050 | Center Action Circular Cup Biopsy Forceps, 50 cm |

### Grasping Forceps

|                          |        |                                                 |
|--------------------------|--------|-------------------------------------------------|
| <input type="checkbox"/> | 505650 | Pediatric Grasping Forceps, 35 cm               |
| <input type="checkbox"/> | 506004 | Forward Grasping Forceps Small, 40 cm           |
| <input type="checkbox"/> | 505020 | Center Action Alligator Grasping Forceps, 50 cm |
| <input type="checkbox"/> | 505032 | Center Action Kahler Grasping Forceps, 50 cm    |
| <input type="checkbox"/> | 505404 | Jackson Alligator Grasping Forceps, 50 cm       |
| <input type="checkbox"/> | 506016 | Forward Grasping Forceps Large, 50 cm           |

### Grasping Forceps

|                          |        |                                                 |
|--------------------------|--------|-------------------------------------------------|
| <input type="checkbox"/> | 505650 | Pediatric Grasping Forceps, 35 cm               |
| <input type="checkbox"/> | 506004 | Forward Grasping Forceps Small, 40 cm           |
| <input type="checkbox"/> | 505020 | Center Action Alligator Grasping Forceps, 50 cm |
| <input type="checkbox"/> | 505032 | Center Action Kahler Grasping Forceps, 50 cm    |
| <input type="checkbox"/> | 505404 | Jackson Alligator Grasping Forceps, 50 cm       |
| <input type="checkbox"/> | 506016 | Forward Grasping Forceps Large, 50 cm           |

### Suctions

|                          |        |                                   |
|--------------------------|--------|-----------------------------------|
| <input type="checkbox"/> | 507023 | Aspirating Tube Open End, 35 cm   |
| <input type="checkbox"/> | 507013 | Aspirating Tube Open End, 45 cm   |
| <input type="checkbox"/> | 507014 | Aspirating Tube Open End, 50 cm   |
| <input type="checkbox"/> | 507031 | Aspirating Tube Velvet Eye, 35 cm |
| <input type="checkbox"/> | 507033 | Aspirating Tube Velvet Eye, 45 cm |
| <input type="checkbox"/> | 507034 | Aspirating Tube Velvet Eye, 50 cm |

## Extra-Long Esophageal Instruments (60 cm)

For use with 50cm Esophagoscopes.

|                          |        |                                           |
|--------------------------|--------|-------------------------------------------|
| <input type="checkbox"/> | 505134 | Roberts Straight 4 mm Cup Forceps, 60 cm  |
| <input type="checkbox"/> | 505146 | Roberts Angled 4 mm Cup Forceps, 60 cm    |
| <input type="checkbox"/> | 505406 | Jackson Alligator Grasping Forceps, 60 cm |
| <input type="checkbox"/> | 505515 | Kahler Grasping Forceps, 60 cm            |
| <input type="checkbox"/> | 507015 | Aspirating Tube Open End, 60 cm           |
| <input type="checkbox"/> | 507035 | Aspirating Tube Velvet Eye, 60 cm         |

## Tray Setup

### Standard-Length Instrument Tray

|                          |           |                                                     |
|--------------------------|-----------|-----------------------------------------------------|
| <input type="checkbox"/> | P37304    | Single-Tier ENT Tray w/o Brackets, 11" x 27" x 3.5" |
| <input type="checkbox"/> | P37306 x2 | Universal Keyhole Bracket 15 Position, 10"          |
| <input type="checkbox"/> | 087396    | Ultrasoft Silicone Spiked Mat, 11" W x 22" L        |


### Extra Long Instrument Tray

For storage of 60 cm Esophageal Instruments.

|                          |        |                                             |
|--------------------------|--------|---------------------------------------------|
| <input type="checkbox"/> | P37319 | Extra-Long Tray w/Brackets, 6.5" x 31" x 3" |
|--------------------------|--------|---------------------------------------------|



# Panendoscopy Checklist



Panendoscopy Configuration in Tray P37304 with Tray Shelf P37338

## Laryngoscopes

Light carrier included.

|                          |        |                                                                    |
|--------------------------|--------|--------------------------------------------------------------------|
| <input type="checkbox"/> | 522190 | Ossoff-Pilling Laser Laryngoscope, 19 mm x 168 mm                  |
| <input type="checkbox"/> | 522221 | Low-Profile Pilling™ Dedo-Style Laser Laryngoscope, 19 mm x 178 mm |

## Bronchoscopes

### Holinger-Type Laser Bronchoscopes

Light carrier and window plug included.

Contain ports for flexible instruments and suction.

|                          |         |                                           |
|--------------------------|---------|-------------------------------------------|
| <input type="checkbox"/> | 523140P | Pilling™ Laser Bronchoscope, 7 mm x 40 cm |
| <input type="checkbox"/> | 523160P | Pilling™ Laser Bronchoscope, 9 mm x 40 cm |

## Esophagoscopes

### Adult Esophagoscopes – Oval

Light carrier included.

|                          |        |                                                         |
|--------------------------|--------|---------------------------------------------------------|
| <input type="checkbox"/> | 523648 | Pillsbury Jesberg Esophagoscope, 10 mm x 14 mm x 29 cm  |
| <input type="checkbox"/> | 523642 | Pillsbury Jesberg Esophagoscope, 6.5 mm x 10 mm x 35 cm |

## Instruments

### Biopsy Cup Forceps

|                          |        |                                             |
|--------------------------|--------|---------------------------------------------|
| <input type="checkbox"/> | 505690 | Flexible Circular Cup Forceps, 2 mm x 50 cm |
| <input type="checkbox"/> | 505100 | Jackson 4 mm Straight, 28 cm                |
| <input type="checkbox"/> | 505104 | Jackson 4 mm Angled Up, 28 cm               |
| <input type="checkbox"/> | 505130 | Roberts 4 mm Straight, 40 cm                |
| <input type="checkbox"/> | 505142 | Roberts 4 mm Angled, 40 cm                  |
| <input type="checkbox"/> | 505132 | Roberts Straight 4 mm Cup Forceps, 50 cm    |
| <input type="checkbox"/> | 505144 | Roberts Angled 4 mm Cup Forceps, 50 cm      |

### Grasping Forceps

|                          |        |                                 |
|--------------------------|--------|---------------------------------|
| <input type="checkbox"/> | 509305 | Ossoff Straight Grasping, 28 cm |
| <input type="checkbox"/> | 505422 | Jackson Delicate, 40 cm         |
| <input type="checkbox"/> | 505414 | Double-Action, 50 cm            |

### Suctions

|                          |        |                                    |
|--------------------------|--------|------------------------------------|
| <input type="checkbox"/> | 507010 | 12 Fr Jackson Open End, 20 cm      |
| <input type="checkbox"/> | 507030 | 12 Fr Jackson Velvet Eye, 20 cm    |
| <input type="checkbox"/> | 507012 | Aspirating Tube, Open End, 40 cm   |
| <input type="checkbox"/> | 507032 | Aspirating Tube, Velvet Eye, 40 cm |
| <input type="checkbox"/> | 507014 | Aspirating Tube Open End, 50 cm    |
| <input type="checkbox"/> | 507034 | Aspirating Tube Velvet Eye, 50 cm  |

## Accessories

### Teeth Guards

|                          |        |                                     |
|--------------------------|--------|-------------------------------------|
| <input type="checkbox"/> | 502275 | Disposable Teeth Guards 24 Pair/Box |
|--------------------------|--------|-------------------------------------|

### Light Cables

|                          |         |                                                   |
|--------------------------|---------|---------------------------------------------------|
| <input type="checkbox"/> | P358PRA | Pilling™ to Pilling™ Blue Fiber Optic Cable, 8 ft |
| <input type="checkbox"/> | A358PRA | ACMI to Pilling™ Blue Fiber Optic Cable, 8 ft     |

### Visualization

|                          |        |                              |
|--------------------------|--------|------------------------------|
| <input type="checkbox"/> | 524050 | Telescope 0°, 4.0 mm X 45 cm |
| <input type="checkbox"/> | 524037 | 4.0 mm Telescope Guide       |
| <input type="checkbox"/> | 502540 | Fluvog Slide Assembly        |


## Tray Setup

|                          |           |                                                     |
|--------------------------|-----------|-----------------------------------------------------|
| <input type="checkbox"/> | P37304    | Single-Tier ENT Tray w/o Brackets, 11" x 27" x 3.5" |
| <input type="checkbox"/> | P37338    | Tray Shelf, 10.5" x 21.5" x 2.25"                   |
| <input type="checkbox"/> | P37545 x4 | "L" Brackets to Support Shelf                       |
| <input type="checkbox"/> | P37306 x5 | Universal Keyhole Bracket 15 Position, 10"          |
| <input type="checkbox"/> | P37379 x2 | Keyhole Bracket, Single, 1.5"                       |
| <input type="checkbox"/> | P37305    | Accessory Box w/Hardware, 1.5" x 3.5" x 1"          |
| <input type="checkbox"/> | 087396 x2 | Ultrasoft Silicone Spiked Mat, 11" W x 22" L        |

# Additional Accessories

## Laryngoscope Holders / Table Mounted Tray

|                          |        |                                                 |
|--------------------------|--------|-------------------------------------------------|
| <input type="checkbox"/> | 502245 | Lewy-style Laryngoscope Holder                  |
| <input type="checkbox"/> | 502248 | Tobey Laryngoscope Holder with Quick Release    |
| <input type="checkbox"/> | 502260 | Boston University Suspension System             |
| <input type="checkbox"/> | 502261 | Boston University Laryngoscope Handle Extension |
| <input type="checkbox"/> | 502265 | Dedo Table Mounted Mayo Tray                    |



Tobey Laryngoscope Holder with Quick Release 502248

# Lightsource and Accessories

## Headlight

| ITEM NUMBER | PRODUCT DESCRIPTION                                          |
|-------------|--------------------------------------------------------------|
| 528339PS    | Variable 20-110 mm Spot with Bifurcated 9' Fused Fiber Cable |

## Lightsources

| ITEM NUMBER | PRODUCT DESCRIPTION                               |
|-------------|---------------------------------------------------|
| 528300PS    | Xenon 300 Watt Digital Lightsource, 5 Port Turret |

## Floorstand

| ITEM NUMBER | PRODUCT DESCRIPTION                        |
|-------------|--------------------------------------------|
| 528304PS    | Floorstand, Including Base/Pole/Wheel Base |



300 Watt Digital Xenon Lightsource

# Fiberoptic Cable Compatibility Guide

| Lightsource Endfitting              | Instrument Endfittings  |                                  |           |                                          |           |                            |
|-------------------------------------|-------------------------|----------------------------------|-----------|------------------------------------------|-----------|----------------------------|
|                                     | ACMI® FEMALE INSTRUMENT | BIFURCATED (2) PILLING™ ROTATING |           | PILLING™ ROTATING BAYONET (PILLING BLUE) |           | STORZ® THREADED INSTRUMENT |
| DESCRIPTION                         | 8' LENGTH               | 6' LENGTH                        | 8' LENGTH | 6' LENGTH                                | 8' LENGTH | 8' LENGTH                  |
| Pilling™ Lightsource Endfitting (P) |                         |                                  |           | P356PRA                                  | P358PRA   |                            |
| ACMI® Lightsource Endfitting (A)    | A358AFA                 |                                  |           |                                          | A358PRA   | A358STA<br>A508STA*        |
| Olympus® Lightsource Endfitting (O) |                         |                                  |           |                                          | O358PRA   |                            |
| Storz® Lightsource Endfitting (S)   |                         |                                  |           |                                          | S358PRA   |                            |
| Wolf® Lightsource Endfitting (W)    |                         | W356PXA                          | W358PXA   |                                          | W358PRA   |                            |

Standard Bundle Size is 3.5 mm.

\* 5.0 mm Bundle Size

## Fiberoptic Cable

Step 1: Select lightsource endfitting

Step 2: Select bundle size

Step 3: Select length

Step 4: Select instrument endfitting



To view Lightsource products go to:  
**myTeleflex.com**



## Additional Needs

[illegible]

## Additional Needs

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

## Recommended Configuration(s) Selected and Approved by:

---

Print Name \_\_\_\_\_

Date \_\_\_\_\_

Signature \_\_\_\_\_

Sales Rep Name

Sales Rep Phone Number

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular access, interventional cardiology and radiology, anesthesia, emergency medicine, surgical, and urology. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit [teleflex.com](https://teleflex.com).

Teleflex is the home of Arrow™, Deknatel™, LMA™, Pilling™, QuikClot™, Rüsch™, UroLift™, and Weck™ – trusted brands united by a common sense of purpose.

### Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

### Regional Offices

**United States:** Phone +1 919 544 8000, Toll Free 866 246 6990, [cs@teleflex.com](mailto:cs@teleflex.com), 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

**Latin America:** Phone +1 919 433 4999, [la.cs@teleflex.com](mailto:la.cs@teleflex.com), 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

**International:** Phone +353 (0)9 06 46 08 00, [orders.intl@teleflex.com](mailto:orders.intl@teleflex.com), Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

**Australia** 1300 360 226

**Austria** +43 (0)1 402 47 72

**Belgium** +32 (0)2 333 24 60

**Canada** +1 (0)800 387 9699

**China** (Shanghai) +86 (0)21 6163 0965

**China** (Beijing) +86 (0)10 6418 5699

**Czech Republic** +420 (0)495 759 111

**France** +33 (0)5 62 18 79 40

**Germany** +49 (0)7151 406 0

**Greece** +30 210 67 77 717

**India** +91 (0)44 2836 5040

**Italy** ++39 0362 58911

**Japan** +81 (0)3 6632 3600

**Korea** +82 2 536 7550

**Mexico** +52 55 5002 3500

**Netherlands** +31 (0)88 00 215 00

**New Zealand** 0800 601 100

**Portugal** +351 22 541 90 85

**Singapore** (SEA non-direct sales countries) +65 6439 3000

**Slovak Republic** +421 (0)3377 254 28

**South Africa** +27 (0)11 807 4887

**Spain** +34 918 300 451

**Switzerland** +41 (0)31 818 40 90

**United Kingdom** +44 (0)1494 53 27 61

**For more information, please visit [teleflex.com](https://teleflex.com).**

Teleflex, the Teleflex logo, Arrow, Deknatel, LMA, Pilling, QuikClot, Rüsch, UroLift, and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates in the U.S. and/or other countries.

Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 04/2024.

© 2024 Teleflex Incorporated. All rights reserved.

MC-008918 Rev 1