



RÜSCH®

Care at home

RÜSCH® FLOCATH QUICK™ HYDROPHILIC INTERMITTENT CATHETERS

Your solution for quick, easy intermittent catheterization

Teleflex is a leading global provider of speciality medical devices for a range of procedures in critical care, home care, and surgery. Our mission is to provide solutions that enable healthcare providers to improve outcomes and enhance patient and provider safety.

The FloCath Quick™ Hydrophilic Intermittent Catheter is designed to help reduce insertion friction and comes in an easy-to-use package.

THE "ALL-IN-ONE" INTERMITTENT CATHETER PROVIDES THE ULTIMATE CONVENIENCE

- Integrated compact, discrete package contains a FloCath Quick™ Hydrophilic Catheter with patented handling sheath and a packet of sterile 0.9% saline
- Handling sheath provides an easy grip, allowing for touchless insertion
- Sheath can also be pushed back and used for extended reach
- Soft eye technology and a gentle tapered tip for comfortable insertion
- Color-coded funnel to easily identify size
- Not made with natural rubber latex
- Available in Straight, Female, and Coude
- Sold individually or as part of a kit



Teleflex®

FLOCATH QUICK™ HYDROPHILIC INTERMITTENT CATHETER

FLOCATH QUICK™ HYDROPHILIC INTERMITTENT CATHETER, STRAIGHT 16" 30/BOX

| PRODUCT # | SIZE |
|-----------|--------|
| 220400060 | 6 Fr. |
| 220400080 | 8 Fr. |
| 220400100 | 10 Fr. |
| 220400120 | 12 Fr. |
| 220400140 | 14 Fr. |
| 220400160 | 16 Fr. |

FLOCATH QUICK™ HYDROPHILIC INTERMITTENT CATHETER, FEMALE 7" 30/BOX

| PRODUCT # | SIZE |
|-----------|--------|
| 220500060 | 6 Fr. |
| 220500080 | 8 Fr. |
| 220500100 | 10 Fr. |
| 220500120 | 12 Fr. |
| 220500140 | 14 Fr. |
| 220500160 | 16 Fr. |

FLOCATH QUICK™ HYDROPHILIC INTERMITTENT CATHETER, COUDÉ 16" 30/BOX

| PRODUCT # | SIZE |
|-----------|--------|
| 220600100 | 10 Fr. |
| 220600120 | 12 Fr. |
| 220600140 | 14 Fr. |
| 220600160 | 16 Fr. |
| 220600180 | 18 Fr. |

FLOCATH QUICK™ HYDROPHILIC INTERMITTENT CATHETER KIT*

The FloCath Quick™ Kits offer the convenience of a small, discrete system while maintaining sterility at home or on the go.


FLOCATH QUICK™ HYDROPHILIC INTERMITTENT CATHETER, STRAIGHT 16" 50/BOX

| PRODUCT # | SIZE |
|-----------|--------|
| 221400060 | 6 Fr. |
| 221400080 | 8 Fr. |
| 221400100 | 10 Fr. |
| 221400120 | 12 Fr. |
| 221400140 | 14 Fr. |
| 221400160 | 16 Fr. |

FLOCATH QUICK™ HYDROPHILIC INTERMITTENT CATHETER, FEMALE 7" 50/BOX

| PRODUCT # | SIZE |
|-----------|--------|
| 221500060 | 6 Fr. |
| 221500080 | 8 Fr. |
| 221500100 | 10 Fr. |
| 221500120 | 12 Fr. |
| 221500140 | 14 Fr. |
| 221500160 | 16 Fr. |

*Kit includes gloves (not made with natural rubber latex), antiseptic skin prep, underpad, metered urine collection bag, and a refuse bag for discrete disposal following catheterization



FloCath Quick™ Hydrophilic Intermittent Catheter Kit

FLOCATH® HYDROPHILIC INTERMITTENT CATHETER

(DOES NOT CONTAIN SALINE SACHET; MUST ADD OWN WATER SOURCE)

FLOCATH® HYDROPHILIC INTERMITTENT CATHETER, STRAIGHT TIP 30/BOX

| PRODUCT # | SIZE |
|-----------|-------------------------|
| 220800060 | 6 Fr. |
| 220800080 | 8 Fr. |
| 220800100 | 10 Fr. |
| 220800120 | 12 Fr. |
| 220800140 | 14 Fr. |
| 220800160 | 16 Fr. |
| 220900140 | 14 Fr. Female Length 7" |

FLOCATH® HYDROPHILIC INTERMITTENT CATHETER, COUDÉ TIP 30/BOX

| PRODUCT # | SIZE |
|-----------|--------|
| 220850100 | 10 Fr. |
| 220850120 | 12 Fr. |
| 220850140 | 14 Fr. |
| 220850160 | 16 Fr. |

TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990

Teleflex, FloCath, FloCath Quick, and Rusch are trademarks or registered trademarks of Teleflex Incorporated or its affiliates. © 2017 Teleflex Incorporated. All rights reserved. MC-001339