



RÜSCH®

Care at home

RÜSCH® MMG H₂O® HYDROPHILIC INTERMITTENT CATHETER CLOSED SYSTEM KITS AND SINGLES

YOUR LIFE, YOUR HEALTH, YOUR CHOICE.

The Rusch® MMG H₂O® Hydrophilic Intermittent Catheter Closed System helps maximize user comfort and may aid in reducing the risk of urethral trauma. The unique hydrophilic coating, activated by the integrated sterile 0.9% saline pouch, allows the catheter to glide easily through the urethra.

The Rusch® MMG H₂O® Hydrophilic Intermittent Catheter Closed System Kit contains sterile gloves (not made with natural rubber latex), underpad, antiseptic skin prep, gauze and refuse bag for convenient disposal.



Teleflex®



☀ AT A GLANCE

- The soft silicone introducer is designed to protect the catheter from bacteria residing in the first few centimeters of the urethra
- The completely closed system incorporates both a thumb grip on the collection bag and a ring pull on the protective cap, allowing for ease of use with limited dexterity
- The unique hydrophilic coating and sterile 0.9% saline combine for a smooth, lubricious coating for ease of insertion
- Patented grip sleeve makes handling easier and allows for touchless catheterization
- Collection bag holds up to 1300 mL and allows for convenient, discrete disposal (kit only)
- This product is not made with natural rubber latex

ORDERING INFORMATION

RÜSCH® MMG H₂O® HYDROPHILIC INTERMITTENT CATHETER CLOSED SYSTEM SINGLES

ITEM CODE	SIZE	CASE QUANTITY
21096080	8 Fr	100
21096100	10 Fr	100
21096120	12 Fr	100
21096140	14 Fr	100
21096160	16 Fr	100

RÜSCH® MMG H₂O® HYDROPHILIC CLOSED SYSTEM CATHETER KITS*

ITEM CODE	SIZE	CASE QUANTITY
20096060	6 Fr.	100
20096080	8 Fr.	100
20096100	10 Fr.	100
20096120	12 Fr.	100
20096140	14 Fr.	100
20096160	16 Fr.	100

*Kits contain sterile gloves (not made with natural rubber latex), underpad, antiseptic skin prep, gauze, and refuse bag for convenient disposal.

TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Teleflex, MMG H₂O, and Rusch are trademarks or registered trademarks of Teleflex Incorporated or its affiliates in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. Revised: 02/2017.

© 2017 Teleflex Incorporated. All rights reserved. MC-003014