




RÜSCH®

Care at home



## RÜSCH® TIEMANN TIP INTERMITTENT CATHETERS

The RÜSCH® Tiemann Tip Intermittent Catheters are designed with a soft, flexible curved tip to help ease catheterization in cases involving an enlarged prostate, strictures, or difficult catheterizations.\* The Tiemann Tip Intermittent Catheters are made with PVC, with each tip color-coded to indicate the French size. The notch on the catheter funnel is used to show the tip position.

RÜSCH® TIEMANN TIP INTERMITTENT CATHETERS		
PRODUCT NUMBER	SIZE	BOX QTY
221800-000120	12 Fr.	50
221800-000140	14 Fr.	50
221800-000160	16 Fr.	50
221800-000180	18 Fr.	50
221800-000200	20 Fr.	50
221800-000220	22 Fr.	50
221800-000240	24 Fr.	50

### AT A GLANCE

- Two staggered eyes
- Color-coded solid tips to indicate size
- Notched funnel to indicate tip position
- Not made with natural rubber latex
- Full range of sizes from 12-24 Fr.

**TO PLACE AN ORDER, CONTACT CUSTOMER SERVICE AT 1.866.246.6990**

Rx only

\*Catheter-Associated Urinary Tract Infections: Definitions and Significance in the Urologic Patient. American Urological Association White Paper.

Teleflex and RÜSCH are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.

© 2018 Teleflex Incorporated. All rights reserved. MC-000971 Rev 1