



RUSCH® URETERAL CATHETERS

Clinical applications

Ureteral catheters have a number of clinical applications such as maintaining drainage, infusion of contrast media for studies and ureteral dilation. Each Rusch Ureteral Catheter is made of plastic and the Luer fitting adapter is enclosed. They are color coded for identification of left or right and are single-use, sterile and X-Ray opaque.

They are 65 cm (26") and are calibrated in 1 cm and 5 cm increments. Whistle Tip, Round Tip, Olive Tip and Blassucci Spiral Tip Catheters all have 2 staggered eyes; the Braasch Bulb has only one eye, and the Cone Tip is open-ended with no eyes. All Rusch Ureteral Catheters are sold in boxes of 10 each (5 red and 5 green).

ORDERING INFORMATION

BLASSUCCI SPIRAL

PRODUCT NUMBER	DESCRIPTION	SIZE
333604	Blassucci Spiral Ureteral	4 Fr.
333605	Blassucci Spiral Ureteral	5 Fr.
333606	Blassucci Spiral Ureteral	6 Fr.



BLASSUCCI SPIRAL

BRAASCH BULB

PRODUCT NUMBER	DESCRIPTION	SIZE
332608	Braasch Bulb Ureteral	4 Fr. shaft; 8 Fr. bulb
332610	Braasch Bulb Ureteral	5 Fr.; 10 Fr. bulb



BRAASCH BULB

CONE TIP

PRODUCT NUMBER	DESCRIPTION	SIZE
334108	Cone Tip Ureteral	5 Fr. shaft; 8 Fr. cone
334110	Cone Tip Ureteral	5 Fr. shaft; 8 Fr. cone
334112	Cone Tip Ureteral	5 Fr. shaft; 12 Fr. cone



CONE TIP

CONNECTOR

PRODUCT NUMBER	DESCRIPTION	QTY
336500	Made of amber latex. Used to connect the ureteral catheters to standard urinary drainage systems.	12/box



CONNECTOR

OLIVE TIP

PRODUCT NUMBER	DESCRIPTION	SIZE
331603	Olive Tip Ureteral	3 Fr.
331604	Olive Tip Ureteral	4 Fr.
331605	Olive Tip Ureteral	5 Fr.
331606	Olive Tip Ureteral	6 Fr.



OLIVE TIP

ROUND TIP

PRODUCT NUMBER	DESCRIPTION	SIZE
332103	Round Tip Ureteral	3 Fr.
332104	Round Tip Ureteral	4 Fr.
332105	Round Tip Ureteral	5 Fr.
332106	Round Tip Ureteral	6 Fr.



ROUND TIP

WHISTLE TIP

PRODUCT NUMBER	DESCRIPTION	SIZE
331103	Whistle Tip Ureteral	3 Fr.
331104	Whistle Tip Ureteral	4 Fr.
331105	Whistle Tip Ureteral	5 Fr.
331106	Whistle Tip Ureteral	6 Fr.



WHISTLE TIP

TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990

Teleflex and Rusch are registered trademarks of Teleflex Incorporated or its affiliates.
©2012 Teleflex Incorporated. All rights reserved. 2012-1428