



RUSCH® BELLY BAG® URINARY COLLECTION DEVICE



The Rüsch® Belly Bag® Urinary Collection Device from Teleflex offers a solution for those using an indwelling catheter following surgery, ambulatory patients, and those leading a more active lifestyle and who prefer a more discrete urine drainage bag. The Belly Bag is worn around the waist with a woven belt and has a quick release buckle.

BENEFITS

- Conceals more effectively than a leg bag
- Enhances comfort and mobility
- Helps to eliminate changing between bed and leg bags
- Helps prevent inadvertent catheter extraction¹

RUSCH® BELLY BAG® URINARY COLLECTION DEVICE

PRODUCT CODE	DESCRIPTION	BOX QTY
B1000	Standard Belly Bag® Collection Device with twist drain	10
B1000CT	Standard Belly Bag® Collection Device with twist drain and 24 extension drain tube	10
B1000P	Standard Belly Bag® Collection Device with twist drain and sample port	10

AT A GLANCE

- 1000 mL capacity
- Cloth backing for added comfort
- Twist valve drain port
- Not made with natural rubber latex
- Sterile
- Anti-Reflux valve
- 52" belt can be trimmed to fit

TO PLACE AN ORDER, CONTACT CUSTOMER SERVICE AT 1.866.246.6990

1. Munnings L.J., Cawood C. Clinical Study of a New Urine Collection Bag. Urologic Nursing. 2003; 23(4):288-291

Teleflex, the Teleflex logo, Belly Bag and Rusch are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this catalogue may not be available in all countries. Please contact your local representative. All data current at time of printing (xx/2016). Subject to technical changes without further notice.

© 2016 Teleflex Incorporated. All rights reserved. MC-001220