



AUTOFUSER®
Disposable Pain Pumps



AutoFuser®



ERGONOMIC PUMP DESIGN

- Rigid, rounded, low-profile housing
- Graduation marks for infusion monitoring
- Not manufactured with natural rubber latex
- Front fill port and flat pump body designed to improve ease of filling
- Reservoir design significantly reduces pressure needed to fill balloon

IN-LINE FILTER

- 0.2 micron bacteria filter
- Air eliminating filter

PAIN PUMP CONFIGURATIONS

Reservoir size, basal rate, and bolus options to meet your needs.

- Two reservoir sizes: 275 mL and 550 mL
- Three fixed basal rates
- Optional Auto Selector variable rate selector
- Optional parallel bolus

PARALLEL BOLUS (OPTIONAL)

- Provides uninterrupted continuous infusion during bolus usage
- Two choices of bolus dosages
- Three choices of bolus refill time

AUTO SELECTOR® INTERFACE (OPTIONAL)

- Two options for variable rate selector:
 - 0.5–7 mL/hr
 - 1–14 mL/hr

For most patients, postoperative pain is their primary concern before surgery and nearly 80% of patients who experience postoperative pain report experiencing moderate to extreme pain.¹ Continuous peripheral nerve blocks have been proven to offer superior postoperative analgesia and improve patient satisfaction compared with opioids, while resulting in fewer minor complications, including nausea/vomiting, pruritis, and sedation².

When used for continuous peripheral nerve block infusions, AutoFuser disposable pain pumps can help effectively manage postoperative pain. AutoFuser pumps offer a combination of a continuous hands-free infusion with an optional on-demand patient controlled bolus to immediately address breakthrough pain.

AUTOFUSER PAIN PUMP ORDERING INFORMATION

FIXED RATE OPTIONS

ORDER #	VOLUME (mL)	FLOW RATE (mL/hr)	BOLUS (mL)	REFILL TIME (min)
MC0020L-CPNB	275	2	none	N/A
MC0040L-CPNB	275	4	none	N/A
MC0050L-CPNB	275	5	none	N/A
MC0050XL-CPNB	550	5	none	N/A

FIXED RATE WITH BOLUS OPTIONS

ORDER #	VOLUME (mL)	FLOW RATE (mL/hr)	BOLUS (mL)	REFILL TIME (min)
MT2060L-CPNB	275	2	2	60
MT4030L-CPNB	275	4	2	30
MT4060L-CPNB	275	4	2	60
MT5015L-CPNB	275	5	2	15
MT5015XL-CPNB	550	5	2	15
MT5060L-CPNB	275	5	2	60
MT5060XL-CPNB	550	5	2	60
MX5060XL-CPNB	550	5	5	60
MT8015XL-CPNB	550	8	2	15
MT8060XL-CPNB	550	8	2	60

AUTO SELECTOR INTERFACE OPTIONS

ORDER #	VOLUME (mL)	FLOW RATE (mL/hr)	BOLUS (mL)	REFILL TIME (min)
MVAL-CPNB	275	0.5–7	none	N/A
MVBL-CPNB	275	1–14	none	N/A
MVBXL-CPNB	550	1–14	none	N/A

AUTO SELECTOR INTERFACE WITH BOLUS OPTIONS

ORDER #	VOLUME (mL)	FLOW RATE (mL/hr)	BOLUS (mL)	REFILL TIME (min)
MVAT60L-CPNB	275	0.5–7	2	60
MVAT60XL-CPNB	550	0.5–7	2	60
MVAX30XL-CPNB	550	0.5–7	5	30
MVAX60XL-CPNB	550	0.5–7	5	60
MVBT15XL-CPNB	550	1–14	2	15
MVBT30XL-CPNB	550	1–14	2	30
MVBT60XL-CPNB	550	1–14	2	60
MVBX30XL-CPNB	550	1–14	5	30
MVBX60XL-CPNB	550	1–14	5	60

TO ORDER, CALL 1-800-788-7999

To learn more about AutoFuser pain pumps and our line of Arrow® peripheral nerve block products, visit teleflex.com/pnb



REFERENCES

1. Warfield, CA, et al. Acute Pain Management: Programs in U.S. Hospitals and Experiences and Attitudes Among U.S. Adults. *Anesthesiology*. 1995; 83: 1090-4.
2. Richman JM, et al. Does Continuous Peripheral Nerve Block Provide Superior Pain Control to Opioids? A Meta-Analysis. *Anesth Analg*. 2006; 102: 248-57.

Arrow and Teleflex are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.
AutoFuser and Auto Selector are registered trademarks of ACE Medical Co., Ltd.

Teleflex is a global provider of medical products designed to enable healthcare providers to protect against infections and improve patient and provider safety. The company specializes in products and services for vascular access, respiratory, general and regional anesthesia, cardiac care, urology and surgery. Teleflex also provides specialty products for device manufacturers.

© 2013 Teleflex Incorporated. All rights reserved. 2013-1763



TELEFLEX MEDICAL PO Box 12600 Research Triangle Park, NC 27709
Toll Free: 866.246.6990 Phone: 919.544.8000 Intl: 919.433.8088
TELEFLEX.COM