



RÜSCH®

Care at Home


RÜSCH® EASYCATH™ INTERMITTENT CATHETERS

RÜSCH® EASYCATH™ INTERMITTENT CATHETER

The EasyCath™ Intermittent Catheter features smooth, polished Soft-Eye™ Technology and a gently tapered tip for a more comfortable, safer insertion. The clear PVC material eliminates any concerns of latex allergies. A full line of Straight Tip, Coudé Tip and Female lengths are available.

RÜSCH EASYCATH INTERMITTENT CATHETER KITS

The full line of EasyCath Kits provides the convenience of a small, discrete system while maintaining complete sterility at home or on the go. The compact kit includes an EasyCath PVC Catheter in Standard, Coudé or Female styles with gloves (not made with natural rubber latex), lubricant, antiseptic skin prep, underpad, metered urine collection bag and a refuse bag.



AT A GLANCE

■ Soft-Eye™ Technology polished, recessed eyelets for enhanced comfort

■ Full range of sizes from 6-16 Fr.

■ Color-coded funnel to easily identify size

Teleflex®



KIT CONTENTS

- Underpad
- One pair of gloves (not made with natural rubber latex)
- Packet of lubricating jelly
- Antiseptic skin prep
- 1500 mL metered collection bag
- Refuse bag

ORDERING INFORMATION

RÜSCH® EASYCATH™ CATHETER STRAIGHT TIP 50/BOX

PRODUCT #	SIZE	LENGTH	PACKAGING
EC060	6 Fr.	11"	Straight
EC080	8 Fr.	11"	Straight
EC100	10 Fr.	11"	Straight
EC120	12 Fr.	16"	Curved
EC121	12 Fr.	16"	Straight
EC140	14 Fr.	16"	Curved
EC141	14 Fr.	16"	Straight
EC160	16 Fr.	16"	Curved
EC161	16 Fr.	16"	Straight
ECS140-Soft PVC	14 Fr.	16"	Straight

RÜSCH® EASYCATH™ CATHETER COUDÉ TIP 50/BOX

PRODUCT #	SIZE	LENGTH	PACKAGING
EC103	10 Fr.	11"	Straight
EC123	12 Fr.	16"	Curved
EC129	12 Fr.	16"	Straight
EC143	14 Fr.	16"	Curved
EC149	14 Fr.	16"	Straight
EC163	16 Fr.	16"	Curved

RÜSCH® EASYCATH™ CATHETER FEMALE 50/BOX

PRODUCT #	SIZE	LENGTH	PACKAGING
EC065	6 Fr.	7"	
EC085	8 Fr.	7"	
EC105	10 Fr.	7"	
EC125	12 Fr.	7"	
EC144	12 Fr.	7"	No Funnel
EC145	14 Fr.	7"	
EC165	16 Fr.	7"	

RÜSCH® EASYCATH™ CATHETER STRAIGHT TIP KIT 50/BOX

PRODUCT #	SIZE
ECK060	6 Fr.
ECK080	8 Fr.
ECK100	10 Fr.
ECK120	12 Fr.
ECK140	14 Fr.
ECK160	16 Fr.

RÜSCH® EASYCATH™ CATHETER COUDÉ TIP KIT 50/BOX

PRODUCT #	SIZE
ECK103	10 Fr.
ECK123	12 Fr.
ECK143	14 Fr.
ECK163	16 Fr.

RÜSCH® EASYCATH™ CATHETER FEMALE LENGTH KIT 50/BOX

PRODUCT #	SIZE
ECK065	6 Fr.
ECK085	8 Fr.
ECK105	10 Fr.
ECK125	12 Fr.
ECK145	14 Fr.
ECK165	16 Fr.

TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990

Rx only

Teleflex, EasyCath, RUSCH and Soft-Eye are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.
© 2015 Teleflex Incorporated. All rights reserved. MC-000064