



RÜSCH®

Care at home



## RÜSCH® TIEMANN TIP INTERMITTENT CATHETERS

The RÜSCH® Tiemann Tip Intermittent Catheters are designed with a soft, flexible curved tip to help ease catheterization in cases involving an enlarged prostate, strictures, or difficult catheterizations.\* The Tiemann Tip Intermittent Catheters are made with PVC, with each tip color-coded to indicate the French size. The notch on the catheter funnel is used to show the tip position.

RÜSCH® TIEMANN TIP INTERMITTENT CATHETERS		
PRODUCT NUMBER	SIZE	BOX QTY
221800120	12 Fr.	50
221800140	14 Fr.	50
221800160	16 Fr.	50
221800180	18 Fr.	50
221800200	20 Fr.	50
221800220	22 Fr.	50
221800240	24 Fr.	50

### AT A GLANCE

- Two staggered eyes
- Color-coded solid tips to indicate size
- Notched funnel to indicate tip position
- Not made with natural rubber latex
- Full range of sizes from 12-24 Fr.

TO PLACE AN ORDER, CONTACT CUSTOMER SERVICE AT 1.866.246.6990

Rx only

\*Catheter-Associated Urinary Tract Infections: Definitions and Significance in the Urologic Patient. American Urological Association White Paper.

Teleflex and RÜSCH are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.

© 2015 Teleflex Incorporated. All rights reserved. MC-000971