

Arrow®



Endurance™ Extended Dwell Peripheral Catheter System

When a PIV is not enough and a PICC is too much.


When a PIV is not enough and a PICC is too much.

The Arrow® Endurance™ Extended Dwell Peripheral Catheter System from Teleflex is designed to provide continuous and reliable IV therapy for the entire length of stay. With a familiar insertion procedure, the Arrow® Endurance™ Catheter provides clinicians with another option when choosing the right line for their patients.

- **88% first attempt success rate¹**
- **Minimized blood exposure:**
 - Closed catheter system designed to enhance clinician safety
- **Minimized kinking:**
 - Variable stiffness guidewire
 - Angled stabilization wings and flexible catheter nose
- **Extended dwell/short term use²**




Closed catheter system designed for added safety.



Arrow® Endurance™ Extended Dwell Peripheral Catheter System

| DESCRIPTION | | DEVICE | | BASIC KIT | |
|-------------|-------------------------------|-----------|-----------|-------------|-------------|
| | | 6 CM | 8 CM | 6CM | 8 CM |
| Blue | 22 Ga Extended Dwell Catheter | EDC-00622 | EDC-00822 | EDC-00622-B | EDC-00822-B |
| Pink | 20 Ga Extended Dwell Catheter | EDC-00620 | EDC-00820 | EDC-00620-B | EDC-00820-B |
| Green | 18 Ga Extended Dwell Catheter | n/a | EDC-00818 | n/a | EDC-00818-B |



Device



Basic Kit

For custom procedural kits that fit your facility's specific preferences, contact your Teleflex sales representative or call 800.523.8446.

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow®, Deknatel®, Hudson RCI®, LMA®, Pilling®, Rüsçh® and Weck® – trusted brands united by a common sense of purpose.

Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

Regional Offices

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Australia 1300 360 226

Austria +43 (0)1 402 47 72

Belgium +32 (0)2 333 24 60

Canada +1 (0)800 387 9699

China (Shanghai) +86 (0)21 6163 0965

China (Beijing) +86 (0)10 6418 5699

Czech Republic +420 (0)495 759 111

France +33 (0)5 62 18 79 40

Germany +49 (0)7151 406 0

Greece +30 210 67 77 717

India +91 (0)44 2836 5040

Italy ++39 0362 58911

Japan +81 (0)3 6632 3600

Korea +82 2 536 7550

Mexico +52 55 5002 3500

Netherlands +31 (0)88 00 215 00

New Zealand 0800 601 100

Portugal +351 22 541 90 85

Singapore (SEA non-direct sales countries) +65 6439 3000

Slovak Republic +421 (0)3377 254 28

South Africa +27 (0)11 807 4887

Spain +34 918 300 451

Switzerland +41 (0)31 818 40 90

United Kingdom +44 (0)1494 53 27 61

For more information, please visit teleflex.com.

Rx only

Reference

1. Teleflex. Unpublished data. Manuscript in progress. May 2017

2. FDA device regulation 21 CFR 880.5200

Teleflex, the Teleflex logo, Arrow, Deknatel, Endurance, Hudson RCI, LMA, Pilling, Rüsçh, and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All other trademarks are trademarks of their respective owners. Information in this document is not a substitute for the product Instructions for Use. The products in this catalogue may not be available in all countries. Please contact your local representative. Revised: (07/2017).

© 2017 Teleflex Incorporated. All rights reserved. MC-001483 Rev 2